

DNA SEQUENCE-BASED HLA TYPING METHOD

Technical Field

5 The present invention relates to a process for determining genotypes of highly polymorphic systems, such as the major histocompatibility complex of humans, including Class I and Class II HLA genes. Specifically, the method of the present invention involves amplifying 10 the alleles carried by any given individual at a gene locus or loci of interest by polymerase chain reaction with conserved and non-conserved oligonucleotide primers. The polymerase chain reaction products are directly sequenced followed by evaluation of the 15 resulting nucleic acid ladders to determine the genotype of sample nucleic acid.

Background of the Invention

20 The major histocompatibility complex (MHC) includes the human leukocyte antigens (HLA) gene complex which is located on the short arm of human chromosome six. This region encodes cell-surface proteins which regulate the cell-cell interactions of the immune response. The various HLA Class I loci encode the HLA 25 antigens, 44,000 dalton polypeptides which associate with B-2 microglobulin at the cell surface. The Class I molecules are involved in the recognition of target cells by cytotoxic T lymphocytes. HLA Class II loci encode cell surface heterodimers, composed of proteins 30 of 29,000 and 34,000 daltons, respectively. These Class II molecules are also involved in the recognition of target cells by helper T lymphocytes.

35 The HLA-A, HLA-B, and HLA-C loci of the HLA Class I region as well as the HLA-DPB, HLA-DQB, HLA-DQA, HLA-DOB and HLA-DPA loci of the HLA Class II region exhibit an extremely high degree of polymorphism. The WHO nomenclature committee for factors of the HLA system

[Marsh and Bodmer, *Immunogenetics*, 31:131 (1990)]

designated 25 alleles of HLA-A (HLA-A-0101, A-0201, etc.), 32 alleles of HLA-B, and 11 alleles of HLA-C, 43 HLA-DQB alleles, 13 HLA-DOB alleles, 8 HLA-DQA alleles, 5 4 HLA-DPA alleles and 19 HLA-DPB alleles. Since this high degree of polymorphism is thought to relate to the function of the HLA molecules, much effort has gone into determining its molecular basis and the functional 10 implications of its polymorphisms (i.e., in transplantation). With the cloning of certain HLA genes this effort has extended to the DNA level.

The Class II genes of the HLA-D region on the short arm of human chromosome six constitute one of the most polymorphic genetic systems known [Bach, *Immunol. Today*, 6:89 (1985)]. The HLA Class II molecules (DR, DQ and DP) are heterodimeric glycoproteins composed of two non-covalently associated chains (alpha and beta) which serve as restricting elements in nominal antigen presentation in the context of self [Zinkernagel and 20 Doherty, *Nature*, 248:701 (1974)] or as foreign antigens in alloresponses [Bach and Van Rood, *N. Engl. J. Med.*, 225:806 (1976)].

Allelic polymorphism of the HLA-D region encoded specificities can be determined by serological methods for phenotyping, mixed lymphocyte cultures using 25 homozygous typing cells, primed lymphocyte testing, determination of restriction fragment length polymorphisms and, more recently, oligotyping [Bach, *Immunol. Today*, 6:118 (1985); Bidwell (1988); Tiercy et al., *Proc. Natl. Acad. Sci. USA*, 85:1198 (1988)]. Present efforts focus largely on the development of molecular approaches to typing, such as RFLP and oligotyping [Bidwell, *Immunol. Today*, 6:118 (1988); Tiercy et al., *Immunol. Today*, 6:118 (1988); Erlich and Bugawan, in *PCR* 35 *Techniques*, H. A. Erlich, ed., Stockton Press, New York (1989)].

The cloning and sequencing of several HLA-DR, DO, and -DP alleles has revealed that their amino acid polymorphisms are located in hypervariable regions of their N-terminal domains, encoded by the second exon of DRB1, DRB3/4/5, DQA1 and DQB1, DPA1 and DPB1 genes [Marsh and Bodmer, *IMMUNOLOGY* (1990); Todd et al., *NATURE*, 322:559 (1987)]. This information has allowed the design of allele-specific oligonucleotides which can be used in the characterization of the known HLA Class II polymorphisms by means of their hybridization to DNA on a solid support (oligomer typing) or for sequencing [Tacey et al., *IMMUNOLOGY* (1988); Erlich and Bugawan, *IMMUNOLOGY* (1989); Todd et al., *IMMUNOLOGY* (1987); Saitki et al., *SCIENCE*, 220:1350 (1985); Nullis and Talcro, *Methods Enzymol.*, 153:335 (1987); Saitki et al., *Nature*, 324:163 (1986); Scharf et al., *SCIENCE*, 233:1076 (1986); Gyllenstein and Erlich, *Proc. Natl. Acad. Sci. USA*, 85:7652 (1988)]. Oligonucleotide typing, although rapid, requires the use of a rather large number of

oligonucleotides for each locus and cannot detect previously unidentified sequence polymorphisms, likely to exist in non-Caucasian populations; further, the approach may not be easily applicable to and may not be practical for the analysis of Class I polymorphisms.

25 Direct sequencing of single-stranded DNA generated by PCR using allele-specific oligonucleotides has been successfully used to examine polymorphism at DQA1 locus [Gyllenstein and Erlich, *IMMUNOLOGY* (1988)]. Application of this approach to DRB genes is, however, problematic due to the strong sequence homology among DRB1, DRB3, DRB4 and DRB5 genes and the presence of up to four different versions of each of these genes in most individuals (isotypic complexity). The very complex ladders generated by direct sequencing make this present process impractical for accurate and rapid determination of HLA types. Thus, direct sequencing of HLA-PCR products has been limited to previous knowledge of the HLA types

carried by a given individual and as such is not suitable for routine HLA typing [Sach, *IMMUNOLOGY* (1985); Bidwell, *IMMUNOLOGY* (1988); Tacey et al., *IMMUNOLOGY* (1988); Erlich and Bugawan, *IMMUNOLOGY* (1989)].

5 Currently, HLA typing is routinely done in connection with many medical procedures, e.g., organ transplantation. Rejection of organ grafts is believed to be diminished if the HLA alleles of donor and recipient are identical. The numerous alleles of HLA genes in the population also make HLA typing useful for paternity testing. However, the currently available techniques are incapable of differentiating among all of the polymorphisms associated with the alleles at Class I and Class II HLA loci. Other drawbacks to current HLA typing are the availability of standard sera necessary to conduct serological tests, the speed of obtaining test results (i.e., XLC takes 5-7 days), and that only the already known HLA types, but not new polymorphisms, are detected by these techniques. In the case of tissue typing in organ transplants and in relatively high volume genetic evaluations, such as paternity testing, the length of time associated with current HLA typing techniques causes unnecessary delay and the results may not be highly accurate.

25 Accordingly, there is a need for a method to determine genomic information in highly polymorphic systems, such as the HLA gene complex, that addresses the limitations imposed by previous methods. That is, in the case of the HLA gene complex, a system that is capable of determining the nucleotide sequences of the genes carried by any given individual without the need to have previous knowledge of his or her HLA types as defined by other methods. Furthermore, the invention avoids the use of oligonucleotides specific for each known allele. The technique we present is rapid, requires the use of only a small number of oligonucleotide primers, and can readily detect new

5

sequence variants unidentifiable with more conventional approaches. This system is exemplified by its applicability to the analysis of Class II as well as Class I and Class III genes and is automatable.

5

Summary of the Invention

The present invention relates to a method for determining the nucleic acid sequence of one or more polymorphic genes of a subject by amplifying and direct sequencing genomic or complementary DNA molecules for each allele at each gene locus to be sequenced using conserved and non-conserved (non-allele-specific) oligonucleotide primers. In a broad sense, the method of the present invention involves sequence-based typing (SST) which provides for unequivocal determination of genetic polymorphism at any genetic locus of interest by direct, simultaneous, sequence analysis of both genomic DNA or expressed (RNA) copies of such a locus. SST can be employed to determine genetic polymorphism at one or more genetic loci of interest, regardless of the complexity of the polymorphism at these loci, including, for example: (1) simple homozygosity or heterozygosity of a unique locus, as exemplified by DQA or the like; (2) isotypic complexity due to multiple, closely related and isoenzymically linked copies of a locus, as exemplified by DRA or the like; and (3) intra-allelic complexity at a locus compounded by interlocus complexity, such as Class I genes or the like. Most known human genetic polymorphisms are of the first, and simplest, type.

Use of the SST method provides overlapping sequence data comprised of only the copies of the locus of interest as is exemplified by each of the types of HLA loci. The SST strategy is designed to ensure selection of a given locus with equal representation of each copy of that locus by equal amplification and direct sequencing of mixtures of both alleles of that locus and direct interpretation of the overlapping

6

sequencing patterns generated by this approach. Thus, providing a method for determining genetic polymorphism at one or more genetic loci of interest which can be employed, for example, in HLA typing, detection,

5 evaluation, and/or characterization of genetic diseases such as, for example, sickle cell anemia, cystic fibrosis, Thalassemia, and the like, and detection, evaluation, and/or characterization of polymorphism in genetic loci associated with various cancers such as p53, Ras, myc, associated with carcinomas, leukemias, sarcomas or the like.

Use of the method according to the present

invention is exemplified by a system providing for rapid and accurate determination of a major histocompatibility complex class genotype of a subject in a sample (e.g., Class I or Class II). Most particularly, the method is directed to determining at least one HLA Class II gene locus including DRB1, DRB3, DRB4, DOB1, DQA1, DPB1 and DPB1 genes. In the case of Class I genotypes, the method is envisioned as being useful to determine A, B, and C loci genes.

To determine a gene locus nucleic acid sequence, polymorphism with the method of the present invention, nucleic acid (RNA or DNA) from a sample is isolated. In the case of RNA, cDNA molecules for each allele of at least one gene locus to be sequenced are synthesized by employing a locus-specific oligonucleotide primer that anneals to a conserved region of each allele of each gene locus. According to the present invention, the sample nucleic acid sequence is determined by amplifying the cDNA molecules or genomic DNA by polymerase chain reaction to generate sufficient product for each allele of each gene locus to be sequenced, with all of the alleles for each gene locus and chromosome to be sequenced being amplified with at least one conserved oligonucleotide primer pair, and at least one of the alleles of each gene locus and chromosome to be

sequenced being amplified with at least one non-conserved oligonucleotide primer and at least one conserved primer; preparing the products of each PCR for sequencing (clean); sequencing directly the products of each polymerase chain reaction product to detect each allele at each gene locus of each chromosome, with an enzyme appropriate for DNA sequencing, such as Taq polymerase and a conserved primer specific for each locus that is sequenced; and analyzing each sequenced product for each locus and primer combination(s) to determine the genotype of the subject.

In a preferred embodiment of the present invention the sequence of each polymerase chain reaction product for each allele of each gene locus is determined by analyzing each nucleic acid single and/or overlapping ladder generated for each directly sequenced polymerase chain reaction product. The analysis is conducted by comparing the nucleotide sequence of each allele of each gene locus sequence to known sequences for each locus, followed by comparing the sequence of each gene locus amplified with the non-conserved/conserved oligonucleotide primer pair to the nucleotide sequence of each allele of the gene locus amplified with a conserved oligonucleotide primer pair. Comparison of nucleic acid ladders for sequenced alleles can be conducted visually or using computer software.

In a preferred embodiment, the process of the invention is automated for use in rapid genotype determinations, including diagnosis of genetic diseases. Automation of the process includes isolating the sample nucleic acid with an RNA/DNA extractor; amplifying the synthesized cDNA molecule or the isolated DNA molecule by polymerase chain reaction using a thermocycler to generate the polymerase chain reaction products; sequencing the polymerase chain reaction products in an automated sequencing apparatus; and analyzing each sequenced polymerase chain reaction product with the

computer having a database with allelic sequence information and the capacity to conduct the appropriate subtraction algorithm for comparing the polymerase chain reaction product sequence for each allele amplified with a conserved oligonucleotide primer pair to the nucleic acid sequence of each allele sequenced with a non-conserved/conserved oligonucleotide primer pair.

The invention further relates to specific groups of oligonucleotide primers useful in the steps of cDNA synthesis, cDNA/genomic DNA amplification by polymerase chain reaction and direct sequencing of the polymerase chain reaction products to determine the nucleotide sequence of each of the alleles at each locus of each chromosome that is amplified. Useful single strand DNA oligonucleotide primers are described in Table 1 herein.

Brief Description of the Drawings

Figure 1A shows a schematic of the cDNA/PCR/Sequencing experiments for DB3 (DB31, DB33, DB34 and DB35), DQ41, DQ41L, D2A1 and DP41 genes.

Figure 1B shows a schematic of the primer binding sites on DB3, DQ41, DQ41L, D2A1 and DP41 transcripts. Stippled boxes represent primers used in the cDNA synthesis reactions; black boxes represent conserved (or Type 1) primers, used for PCR; checked boxed represent non-conserved (or Type 2) primers, also used for PCR; and blank boxes represent sequencing primers.

Figure 1C shows a schematic of the primer binding sites on DQ41, DB41, DP41 and DP41 genes in their germline configuration. Only these primers exclusively used for genomic DNA samples are shown in the Figure.

Figure 2A shows a flow-chart of the procedure for peripheral blood samples. Each reaction is performed in a different test tube. The reactions are

9

named with capital letter in parenthesis; these letters correspond to those shown in Table II (combinations of primers/reaction). Only the "routine" combinations of primers are shown in this Figure.

Figure 2B is a flow-chart of the procedure for forensic samples, where DNA is usually the only available genetic material to work with. DNA in these situations is usually isolated from hair, sperm, blood stains, etc. The combinations of primers per reaction shown in the Figure correspond to the "routine" combinations only.

Figure 3 shows direct sequencing of Class II HLA d2DNA generated using conserved oligonucleotides.

Lanes are read from left to right as G-A-T-C. 1, DQB1 ladder for a DQB1*0201/DQB1*0302 heterozygote; 2, DQB1 ladder for a DQB1*0103 homozygous cell line; 3, DRB1 ladder for a DRB1*0301, DRB3*0101/DRB1*0401, DRB4*0101 heterozygote. Positions where there is more than one band are indicated on the side of the ladder and the templates they correspond to are indicated at the top of the Figure. To read unambiguously the last 50-60 base pairs of the ladder it is necessary to electrophorese the sequencing gel for an additional hour. Note that the ladders corresponding to the genes at DRB3 or DRB4 loci are fainter in comparison to those corresponding to the genes at DRB1 locus, possibly due to their lower levels of expression. These differences in intensity are generally reproducible and help read the complex overlapping patterns. The positions of the first exon base pair and codon (in parenthesis) that can be read in this Figure are indicated at the bottom of each ladder.

Figure 4 shows direct sequencing of Class II HLA DRB1 cDNA generated using non-conserved oligonucleotides. Lanes are read from left to right as G-A-T-C. Lane 1, DRB1*0101/DRB1*1501, DRB5*0101 heterozygote cDNA amplified with primer DRB17 (selects DRB5*0101 cDNA) (left) and DRB16 (selects DRB1*0101

10

cDNA) (right); lane 2, DRB1*1405, DRB4*0101/DRB1*0301, DRB3*0101 heterozygote cDNAs amplified with the 5' primers DRB17 (selects DRB1*0301 and DRB3*0101 cDNAs) (left) and DRB16 (selects DRB1*0405 and DRB4*0101) (right). Positions where there is more than one band or where the two ladders generated with each primer differ are indicated on the side of the ladders.

Detailed Description of the Invention

As used herein, the term "gene" refers to a segment of DNA, composed of a transcribed region and a regulatory sequence that makes possible a transcription. The term "gene locus" refers to the specific place on the chromosome where a gene is located. The term "allele" refers to the multiple forms of a gene that can exist at a single gene locus at a single chromosome and are distinguishable from the other possible alleles by their differing effects on phenotype (detectable outward manifestations of a specific genotype). "Haplotype" refers to the specific allele composition of the genes at multiple loci on the same chromosome. As used herein the term "genotype" refers to the specific allelic composition of a gene at multiple linked loci at each chromosome (2 haplotypes).

The term "oligonucleotide" as used herein refers to a molecule having two or more deoxyribonucleotides or ribonucleotides, preferably more than three deoxyribonucleotides. The exact number of nucleotides in the molecule will depend on the function of the specific oligonucleotide molecule. As used herein the term "primer" refers to a single stranded DNA oligonucleotide sequence, preferably produced synthetically which is capable of acting as a point of initiation for synthesis of a primer extension product which is complementary to a nucleic acid strand to be copied or a point of initiation for sequencing a DNA molecule. In the case of primers intended for use in

11

synthesizing cDNA or amplifying cDNA or genomic DNA molecules by polymerase chain reaction products, the length and sequence of the primer must be sufficient to prime the synthesis of extension products in the presence of a polymerization enzyme. Preferably, the length of the primer is from about 5-50 nucleotides, more preferably from about 5-20 nucleotides. Specific length and sequence of the primer will depend on complexity of required DNA or RNA target templates, as well as conditions of primer employment such as temperature, ionic strength, and MgCl₂ concentration. When nested primers are used for sequencing, the number of base pairs separating the amplification and sequencing primers on the DNA template are also important considerations.

As used herein, "conserved oligonucleotide primer" (Type 1) refers to an oligonucleotide molecule that corresponds to a region of high DNA sequence conservation (i.e., less than 1-2 nucleotide variations). While the conserved primer need not correspond exactly to the nucleotide template to which it anneals, the conserved primer will have minimal, preferably less than one mismatch with the target nucleotide template. Functionally, conserved primers are capable of equally priming the target nucleotide (cDNA, PCR product, etc.) at high stringency conditions. In contrast to this, as used herein, "non-conserved oligonucleotide primer" (Type 2) refers to an oligonucleotide molecule that has an intended number of mismatches with the possible target nucleotide sequences. The intended number of mismatches can vary with a preferred number of mismatches being about 1-12. Non-conserved primers are characterized by their selective binding to a limited number of alleles at a given locus or at a group of highly homologous loci. The non-conserved primer will bind to the more complementary allele or group of alleles (two or less than two) (i.e., fewer number of

12

mismatches between primer and target template sequence). The specific combinations of conserved and non-conserved primers and the number of reactions per locus or loci used herein are specifically designed to obtain highly accurate results with minimal expenditure of time and cost.

The present invention is directed to a process for determining the sequences of the alleles of polymorphic gene systems carried by any given individual, such as, for example, the human HLA system, genes related to different human genetic disorders, such as sickle cell anemia, cystic fibrosis, or the like, as well as gene systems associated with various cancers, such as p53, myc, or the like. The present invention is exemplified by its utility for determining polymorphism at HLA loci, particularly Class II and Class I genes, the most polymorphic human genetic loci known today, using enzymatic amplification and direct sequencing of the gene cDNA molecules using a limited number of primers and avoiding the use of allele specific oligonucleotides as much as possible. The present method is particularly well suited to determining allelic sequences of Class II HLA genes, thereby providing complete HLA Class II genotype information for a subject. Using the method of the present invention complete Class II HLA typing (DR, DQ and DP) can be performed in about 16 to 24 hours or less.

Generally, the method of the present invention involves: extraction of sample nucleic acid; in the case of RNA, generation of cDNA; cDNA or genomic DNA amplification; direct sequencing of amplification products; and analysis of the direct sequence information. Generation of cDNA, amplifying the cDNA and direct sequencing the cDNA amplification products is accomplished using oligonucleotide primers with specific characteristics, such as those described herein.

A. Oligonucleotide Primers

The oligonucleotide primers of the present invention can be synthesized using any known suitable method, such as phosphotriester and phosphodiester methods. Narang et al., *Methods Enzymol.*, 58:190 (1979); Brown et al., *Methods Enzymol.*, 58:109 (1979).

Oligonucleotides can be prepared using a modified solid support such as a Biossearch 8750 DNA synthesizer.

Useful primers can also be isolated from a biological source using appropriate restriction endonucleases which cut double stranded DNA at or near a nucleotide sequence of interest for use as a primer.

B. Extraction of Sample Nucleic Acid

In the process of the present invention any source of nucleic acid can be used as the sample nucleic acid, as long as the sample contains the nucleic acid sequence of interest. For example, the sample chosen for the present method can be RNA, DNA or a DNA/RNA hybrid. Typical samples include peripheral blood mononuclear cells, (PBMC's), lymphoblastoid cell lines (LCL's), hair cells or the like. For determining human HLA Class II and Class I gene polymorphisms LCL's or PBMC's are preferred. The nucleic acid to be isolated

(e.g. RNA or DNA) will depend on the source of genetic material (blood stain, hair, or peripheral blood cells). However, in the case of HLA Class II genes including DRB1-5, DQB1, DQA1, DPB1 the preferred isolated nucleic acid is total cellular RNA when the typing is to be done for transplantation purposes or paternity testing.

For forensic uses, genomic DNA may be the preferred genetic material in which case different primer considerations would be used. Cytoplasmic and poly(A) + RNA can also be used. It is envisioned that isolation of sample nucleic acid for the present process can be automated using a DNA/RNA extractor (such as

Model 341 DNA extractor available from Applied Biosystems, Inc.; Foster City, CA).

C. Generation of cDNA

Complementary DNA (cDNA) of the sample nucleic acid is generated using specific oligonucleotide primers and cloned reverse transcriptase following general conditions suggested by the enzyme manufacturer (Bethesda Research Laboratories, Gaithersburg, MD). Specific differences in type and amount of primers used, dNTP concentrations and elongation times will be readily apparent to those of skill in the art based on the examples that follow.

D. Polymerase Chain Reaction

Amplification of cDNA or genomic DNA for each gene locus of interest is accomplished using the polymerase chain reaction (PCR) as generally described in U.S. Pat. Nos. 4,683,195 and 4,683,202 to Mullis. The PCR consists of many repetitions of a cycle which consists of: (a) a denaturation step, which melts both strands of a DNA molecule; (b) an annealing step, which is aimed at allowing the primers to anneal specifically to the melted strands of the DNA molecule; and (c) an extension step, which incorporates to the primers deoxyribonucleotides complementary to those of the strand of DNA to which the primers are annealed. The PCR process, as indicated in the Examples, can be conducted using a Thermocycler (Perkin-Elmer, Cetus, Emeryville, CA).

The present invention introduces the use of non-conserved oligonucleotides in the PCR procedure specifically designed to solve the problems associated with, for example, detecting, evaluating, and/or characterizing polymorphism at a polymorphic gene locus or loci of an individual. In the case of HLA typing, the use of non-conserved oligonucleotides addresses the

15

problems one would face in performing HLA typing by sequencing DNA amplified exclusively by using conserved oligonucleotides or allele specific oligonucleotides (see below).

- 5 It is understood that the PCR process is designed for the amplification of specific genes with the use of oligonucleotides specific for the particular gene to be amplified. However, even using completely matched primers, in most cases the PCR is not absolutely specific. In the case of HLA typing, for HLA-DNS genes and Class I genes, the use of conserved primers in PCR will generate complex mixtures of templates, which upon direct sequencing will be seen as overlapping sequencing ladders, cumbersome to interpret. Therefore, genes for which the exact nucleotide sequence information is unknown can not be achieved with an adequate level of certainty. Use of non-conserved oligonucleotides which can selectively anneal under high stringency conditions to two or fewer alleles of a gene locus or group of homologous loci can provide sequence information for the different genes at highly homologous loci in complex heterozygote combinations. Thus, the present invention provides a method useful for determining the genotype for polymorphic gene loci. This is of particular importance to HLA typing, and is applicable to Class I HLA typing as well as Class II typing.

The difference between non-conserved primers and allele-specific oligonucleotides resides in that the latter can only be used when the presence of a particular allele is known, and also requires the use of a specific primer for each of the alleles of the polymorphic system. Thus, combining use of a non-conserved primer and conserved primers to amplify the separate alleles of highly homologous polymorphic genes loci can provide simpler DNA polymerase chain reaction product combinations sufficient to allow unambiguous interpretation of direct sequencing ladders of each

16

alleles for genotype determination with moderate expenditure of time and economical cost.

- The conditions used for the PCR reactions are preferably the same except for the temperature used in the annealing step, which is different depending on the type of primer used, conserved (Type 1) or non-conserved (Type 2). Reactions that use the former primer type are preferably performed at 37°C in the annealing step of the cycle, whereas this step is preferably performed at about 55°C to 60°C in reactions that use the latter type of primers. The concentrations of primers, and buffers used will be apparent from and include the process parameters described in the Examples that follow.

15 E. Direct Sequencing of PCR Products

Direct sequencing of double-stranded DNA generated by the PCR is accomplished using an enzyme appropriate to DNA sequencing, such as Taq polymerase, or the like, and specific combinations of reagents at appropriate concentrations. The sequencing procedure can be conducted in an automatic sequencing apparatus such as the 373A Model DNA Sequencer from Applied Biosystems Inc. (Foster City, CA). The reagents, including sequencing primers and nucleic acid termination mixtures will be understood by those of skill in the art based on the direct sequencing procedure specified in the following Examples.

F. Analysis of Direct Sequenced PCR Products

- The nucleic acid ladders resulting from direct sequencing the cDNA or genomic DNA for each gene locus of interest can be assessed visually from autoradiograms or by employing a computer programmed with nucleotide sequence information for all alleles of all haplotypes and procedures for comparing sequenced alleles and known alleles of gene loci of interest. In a preferred embodiment of the present invention, the evaluation of

gene locus alleles involves a two step process: (a) comparison of the gene sequences of each polymerase chain reaction product (i.e., conserved and non-conserved primer products) with a library of known genotype information such as the information obtained on homologous cell lines very well characterized by methods other than sequencing (Marsh and Bodmer, Immunogenetics, 31:131 (1990)) as well as sequences of individual alleles, followed by (b) comparison of direct sequence information for the polymerase chain reaction product of an allele of a gene locus amplified with a conserved oligonucleotide primer pair and polymerase chain reaction product of alleles of a gene locus or loci amplified with a conserved/non-conserved primer pair.

15 This comparison employs a substitution algorithm or visual cancellation of duplicative sequence ladder information to generate the specific sequence information for each allele of a gene locus.

It is envisioned that the process of the present invention can be used to amplify and sequence known and unknown highly polymorphic systems (e.g., genetic disease-related genes, cancer-related genes, and HLA typing, including Class I, Class II, and Class III HLA typing, and the like). The present process is believed to be useful for paternity testing and forensic medicine, with more accuracy than restriction fragment length polymorphism (RFLP), DNA fingerprinting or dot blot-detection systems. While in the latter only a hybridization pattern is observed, direct sequencing of amplified products shows the exact nucleotide sequence of the amplified genes, and hence is more accurate and reliable.

The method is particularly well suited for Class II HLA typing, reducing its costs, increasing its speed and especially improving its accuracy. As evidenced by the following Examples, sequence polymorphism analysis of DRB1, DRB3, DRB4, DRB5, DRB*2,

DRB1, DRB3 and DRB4 genes can be rapidly performed in any subject of unknown HLA type by means of enzymatic amplification and direct sequencing of Class II genes using a limited number of conserved and non-conserved oligonucleotides. The approach described herein is entirely amenable using currently available technology and, as opposed to previously described methods using oligonucleotide probes and dot blots, has the advantage of detecting the presence of new allelic sequences or sequence microheterogeneity at the population level. The methodology of present invention is envisioned to be useful for detailed analyses of the effects of sequence allelism at different Class II HLA loci on graft survival after allogeneic transplantation. The method of the present invention allows rapid and precise sequence analysis of Class II HLA polymorphism in studies of human disease and may be of interest in the search for new Class II sequence variants in large populations of subjects.

20 The present invention is further described by illustration in the following Examples which are not intended to limit the invention.

EXAMPLE I

25 1. Preparation of oligodeoxyribonucleotide Primers and Sequences for Primer Construction for PCR Amplification of Class II HLA Genes

All of the oligodeoxyribonucleotide primers described herewith were synthesized as described below:

Automated Synthesis of oligodeoxyribonucleotide primers: The b-cyanoethylphosphonidites, obtained from Milligen-Bioscience (Norwalk, CT), were sequentially condensed to a nucleoside derivatized controlled pore glass support using a Bioscience 8750 DNA synthesizer. Condensation cycles included deprotection with dichloroacetic acid in dichloromethane, followed by condensation with benzotriazole and capping with acetic

19

anhydride and 3-methylimidazole in tetrahydrofuran and pyridine, with each cycle time being approximately 9 minutes. Yields at each step were >98% as determined by measuring dimethoxytrityl alcohol release. The methodology for oligodeoxynucleotide synthesis is described in Caruthers, et al., *Methods Enzymol.*, 154:287 (1987).

Deprotection and purification of

oligodeoxyribonucleotide primers: Deprotection and purification of oligodeoxyribonucleotide primers was performed using the procedure described by Schulhof et al., *Nucl. Acids Res.*, 15:397 (1987). Briefly, the oligodeoxyribonucleotide was removed from the solid support by exposure to concentrated ammonium hydroxide at room temperature for about one hour. The solution containing the partially deprotected

oligodeoxyribonucleotide was brought to 65°C for 16 hours. Ammonia was removed and the residue was subjected to chromatography on a C18 reverse-phase column (RP 304, BioRad, Richmond, VA) using a linear gradient of 14 to 20% acetonitrile in 0.1 molar

ammonium/triethylamine, pH 7.0. The dimethoxystyryl group was removed from the HPLC-purified oligodeoxyribonucleotide by treatment with 70% acetic acid. The deprotected oligodeoxyribonucleotide was recovered by precipitation in ether, vacuum centrifuged until dry, resuspended in water and quantitated by measuring its absorbance at 260nm.

Using the above procedure, the following oligonucleotide primers corresponding to specified regions of HNA Class II DQA, DQB, DRB, and DPA loci were synthesized (see Table I below) and extensively tested:

TABLE I
Oligonucleotides Used For The cDNA/PCR/Sequencing Reactions

Accession	Gene Type	Length (nt)	Accession
U007	RNA	105-111	D081
U008	RNA	105-111	D081
U009	RNA	105-111	D081
U010	RNA	105-111	D081
U011	RNA	105-111	D081
U012	RNA	105-111	D081
U013	RNA	105-111	D081
U014	RNA	105-111	D081
U015	RNA	105-111	D081
U016	RNA	105-111	D081
U017	RNA	105-111	D081
U018	RNA	105-111	D081
U019	RNA	105-111	D081
U020	RNA	105-111	D081
U021	RNA	105-111	D081
U022	RNA	105-111	D081
U023	RNA	105-111	D081
U024	RNA	105-111	D081
U025	RNA	105-111	D081
U026	RNA	105-111	D081
U027	RNA	105-111	D081
U028	RNA	105-111	D081
U029	RNA	105-111	D081
U030	RNA	105-111	D081
U031	RNA	105-111	D081
U032	RNA	105-111	D081
U033	RNA	105-111	D081
U034	RNA	105-111	D081
U035	RNA	105-111	D081
U036	RNA	105-111	D081
U037	RNA	105-111	D081
U038	RNA	105-111	D081
U039	RNA	105-111	D081
U040	RNA	105-111	D081
U041	RNA	105-111	D081
U042	RNA	105-111	D081
U043	RNA	105-111	D081
U044	RNA	105-111	D081
U045	RNA	105-111	D081
U046	RNA	105-111	D081
U047	RNA	105-111	D081
U048	RNA	105-111	D081
U049	RNA	105-111	D081
U050	RNA	105-111	D081
U051	RNA	105-111	D081
U052	RNA	105-111	D081
U053	RNA	105-111	D081
U054	RNA	105-111	D081
U055	RNA	105-111	D081
U056	RNA	105-111	D081
U057	RNA	105-111	D081
U058	RNA	105-111	D081
U059	RNA	105-111	D081
U060	RNA	105-111	D081
U061	RNA	105-111	D081
U062	RNA	105-111	D081
U063	RNA	105-111	D081
U064	RNA	105-111	D081
U065	RNA	105-111	D081
U066	RNA	105-111	D081
U067	RNA	105-111	D081
U068	RNA	105-111	D081
U069	RNA	105-111	D081
U070	RNA	105-111	D081
U071	RNA	105-111	D081
U072	RNA	105-111	D081
U073	RNA	105-111	D081
U074	RNA	105-111	D081
U075	RNA	105-111	D081
U076	RNA	105-111	D081
U077	RNA	105-111	D081
U078	RNA	105-111	D081
U079	RNA	105-111	D081
U080	RNA	105-111	D081
U081	RNA	105-111	D081
U082	RNA	105-111	D081
U083	RNA	105-111	D081
U084	RNA	105-111	D081
U085	RNA	105-111	D081
U086	RNA	105-111	D081
U087	RNA	105-111	D081
U088	RNA	105-111	D081
U089	RNA	105-111	D081
U090	RNA	105-111	D081
U091	RNA	105-111	D081
U092	RNA	105-111	D081
U093	RNA	105-111	D081
U094	RNA	105-111	D081
U095	RNA	105-111	D081
U096	RNA	105-111	D081
U097	RNA	105-111	D081
U098	RNA	105-111	D081
U099	RNA	105-111	D081
U100	RNA	105-111	D081
U101	RNA	105-111	D081
U102	RNA	105-111	D081
U103	RNA	105-111	D081
U104	RNA	105-111	D081
U105	RNA	105-111	D081

2. Combinations of Primers for cDNA/PCR/Sequence Reactions

There are specific combinations of

- oligonucleotide primers for each reaction and for each locus, including cDNA synthesis, PCR amplification and direct sequencing, which are designed to provide all the necessary sequence information for obtaining highly accurate, fast and inexpensive typing results. These combinations are listed in Table II below as "routine" combinations. In addition, Table II includes a list of "alternative" combinations of oligonucleotides for each locus which may be used to confirm results obtained with the "routine" combinations for a particular locus not expected according to, for instance, known haplotypic maps. These "unexpected" results are usually indicative of the existence of new alleles and/or haplotypes, which can be confirmed with the use of the alternative combinations of oligonucleotides. In any case, each of these combinations of oligonucleotides is characterized by its ability to generate an end-product (sequencing ladder) which is suitable of being accurately read by the naked eye or processed by computer operated under appropriate software.

- For typing purposes in the clinical setting, such as in transplantation, the method uses RNA isolated from peripheral blood mononuclear cells as starting material; for forensic purposes, however, DNA is often the only available template. Although for each template (RNA or DNA) different combinations of oligonucleotides are used (see Table II), the general strategy for typing, including the interpretation of the results is essentially the same. The specific combinations of primers for "routine" RNA and DNA analysis, respectively, are described below in more detail. The general overview of the RNA typing strategy is shown in Figures 1 and 2 and discussed further in Examples 2 and 3.

TABLE II

Combinations of Primers for cDNA/PCR/Seq Reactions

1. RNA	2. DNA	PCR	A.L.	Seq*
Routine	Routine	DB11	37°C	
DB30/DB12/DB42**	DB20	DB11		DB430/DB412
B. 2 DB30	DB20			DB430/DB412
C. 2 DB30	DB20			DB430/DB412
D. 2 DB30	DB20			DB430/DB412
E. 1 DB30	DB20			DB430/DB412
F. 1 DB30	DB20			DB430/DB412
G. 1 DB30	DB20			DB430/DB412
H. 1 DB30	DB20			DB430/DB412
I. 1 DB30	DB20			DB430/DB412
Alternative	Alternative	DB15	55°C	
J. 2 DB30	DB20	DB15		DB430/DB412
K. 2 DB30	DB20	DB15		DB430/DB412
L. 2 DB30	DB20	DB15		DB430/DB412
M. 1 DB30	DB20	DB15		DB430/DB412
N. 1 DB30	DB20	DB15		DB430/DB412
O. 2 DB30	DB20	DB15		DB430/DB412
P. 2 DB30	DB20	DB15		DB430/DB412
Q. 1 DB30	DB20	DB15		DB430/DB412
R. 1 DB30	DB20	DB15		DB430/DB412
2. DNA	2. DNA	PCR2	A.L.	Seq
Routine	Routine	DB22	37°C	
S. 2 DB1405	DB1405	DB22		DB1312/DB1400
T. 2 DB1405	DB1405	DB22		DB1312/DB1400
U. 2 DB1405	DB1405	DB22		DB1312/DB1400
V. 2 DB1405	DB1405	DB22		DB1312/DB1400
W. 1 DB1405	DB1405	DB22		DB1312/DB1400
X. 1 DB1405	DB1405	DB22		DB1312/DB1400
Y. 1 DB1405	DB1405	DB22		DB1312/DB1400
Z. 1 DB1405	DB1405	DB22		DB1312/DB1400
Alternative	Alternative	DB24	37°C	
AA. 2 DB1405	DB1405	DB24		DB1312/DB1400
AB. 2 DB1405	DB1405	DB24		DB1312/DB1400
AC. 2 DB1405	DB1405	DB24		DB1312/DB1400
AD. 1 DB1405	DB1405	DB24		DB1312/DB1400
AE. 1 DB1405	DB1405	DB24		DB1312/DB1400
AF. 1 DB1405	DB1405	DB24		DB1312/DB1400
AG. 1 DB1405	DB1405	DB24		DB1312/DB1400
AH. 1 DB1405	DB1405	DB24		DB1312/DB1400
AI. 1 DB1405	DB1405	DB24		DB1312/DB1400
AJ. 1 DB1405	DB1405	DB24		DB1312/DB1400
AK. 1 DB1405	DB1405	DB24		DB1312/DB1400
AL. 1 DB1405	DB1405	DB24		DB1312/DB1400
AM. 1 DB1405	DB1405	DB24		DB1312/DB1400
AN. 1 DB1405	DB1405	DB24		DB1312/DB1400
AO. 1 DB1405	DB1405	DB24		DB1312/DB1400
AP. 1 DB1405	DB1405	DB24		DB1312/DB1400
AQ. 1 DB1405	DB1405	DB24		DB1312/DB1400
AR. 1 DB1405	DB1405	DB24		DB1312/DB1400
AS. 1 DB1405	DB1405	DB24		DB1312/DB1400
AT. 1 DB1405	DB1405	DB24		DB1312/DB1400
AU. 1 DB1405	DB1405	DB24		DB1312/DB1400
AV. 1 DB1405	DB1405	DB24		DB1312/DB1400
AW. 1 DB1405	DB1405	DB24		DB1312/DB1400
AX. 1 DB1405	DB1405	DB24		DB1312/DB1400
AY. 1 DB1405	DB1405	DB24		DB1312/DB1400
AZ. 1 DB1405	DB1405	DB24		DB1312/DB1400
BA. 1 DB1405	DB1405	DB24		DB1312/DB1400
BB. 1 DB1405	DB1405	DB24		DB1312/DB1400
BC. 1 DB1405	DB1405	DB24		DB1312/DB1400
BD. 1 DB1405	DB1405	DB24		DB1312/DB1400
BE. 1 DB1405	DB1405	DB24		DB1312/DB1400
BF. 1 DB1405	DB1405	DB24		DB1312/DB1400
BG. 1 DB1405	DB1405	DB24		DB1312/DB1400
BH. 1 DB1405	DB1405	DB24		DB1312/DB1400
BI. 1 DB1405	DB1405	DB24		DB1312/DB1400
BJ. 1 DB1405	DB1405	DB24		DB1312/DB1400
BK. 1 DB1405	DB1405	DB24		DB1312/DB1400
BL. 1 DB1405	DB1405	DB24		DB1312/DB1400
BM. 1 DB1405	DB1405	DB24		DB1312/DB1400
BN. 1 DB1405	DB1405	DB24		DB1312/DB1400
BO. 1 DB1405	DB1405	DB24		DB1312/DB1400
BP. 1 DB1405	DB1405	DB24		DB1312/DB1400
BQ. 1 DB1405	DB1405	DB24		DB1312/DB1400
BR. 1 DB1405	DB1405	DB24		DB1312/DB1400
BS. 1 DB1405	DB1405	DB24		DB1312/DB1400
BT. 1 DB1405	DB1405	DB24		DB1312/DB1400
BU. 1 DB1405	DB1405	DB24		DB1312/DB1400
BV. 1 DB1405	DB1405	DB24		DB1312/DB1400
BW. 1 DB1405	DB1405	DB24		DB1312/DB1400
BX. 1 DB1405	DB1405	DB24		DB1312/DB1400
BY. 1 DB1405	DB1405	DB24		DB1312/DB1400
BZ. 1 DB1405	DB1405	DB24		DB1312/DB1400
CA. 1 DB1405	DB1405	DB24		DB1312/DB1400
CB. 1 DB1405	DB1405	DB24		DB1312/DB1400
CC. 1 DB1405	DB1405	DB24		DB1312/DB1400
CD. 1 DB1405	DB1405	DB24		DB1312/DB1400
CE. 1 DB1405	DB1405	DB24		DB1312/DB1400
CF. 1 DB1405	DB1405	DB24		DB1312/DB1400
CG. 1 DB1405	DB1405	DB24		DB1312/DB1400
CH. 1 DB1405	DB1405	DB24		DB1312/DB1400
CI. 1 DB1405	DB1405	DB24		DB1312/DB1400
CJ. 1 DB1405	DB1405	DB24		DB1312/DB1400
CK. 1 DB1405	DB1405	DB24		DB1312/DB1400
CL. 1 DB1405	DB1405	DB24		DB1312/DB1400
CM. 1 DB1405	DB1405	DB24		DB1312/DB1400
CN. 1 DB1405	DB1405	DB24		DB1312/DB1400
CO. 1 DB1405	DB1405	DB24		DB1312/DB1400
CP. 1 DB1405	DB1405	DB24		DB1312/DB1400
CQ. 1 DB1405	DB1405	DB24		DB1312/DB1400
CR. 1 DB1405	DB1405	DB24		DB1312/DB1400
CS. 1 DB1405	DB1405	DB24		DB1312/DB1400
CT. 1 DB1405	DB1405	DB24		DB1312/DB1400
CU. 1 DB1405	DB1405	DB24		DB1312/DB1400
CV. 1 DB1405	DB1405	DB24		DB1312/DB1400
CW. 1 DB1405	DB1405	DB24		DB1312/DB1400
CX. 1 DB1405	DB1405	DB24		DB1312/DB1400
CY. 1 DB1405	DB1405	DB24		DB1312/DB1400
CA. 1 DB1405	DB1405	DB24		DB1312/DB1400
CB. 1 DB1405	DB1405	DB24		DB1312/DB1400
CC. 1 DB1405	DB1405	DB24		DB1312/DB1400
CD. 1 DB1405	DB1405	DB24		DB1312/DB1400
CE. 1 DB1405	DB1405	DB24		DB1312/DB1400
CF. 1 DB1405	DB1405	DB24		DB1312/DB1400
CG. 1 DB1405	DB1405	DB24		DB1312/DB1400
CH. 1 DB1405	DB1405	DB24		DB1312/DB1400
CI. 1 DB1405	DB1405	DB24		DB1312/DB1400
CJ. 1 DB1405	DB1405	DB24		DB1312/DB1400
CK. 1 DB1405	DB1405	DB24		DB1312/DB1400
CL. 1 DB1405	DB1405	DB24		DB1312/DB1400
CM. 1 DB1405	DB1405	DB24		DB1312/DB1400
CN. 1 DB1405	DB1405	DB24		DB1312/DB1400
CO. 1 DB1405	DB1405	DB24		DB1312/DB1400
CP. 1 DB1405	DB1405	DB24		DB1312/DB1400
CQ. 1 DB1405	DB1405	DB24		DB1312/DB1400
CR. 1 DB1405	DB1405	DB24		DB1312/DB1400
CS. 1 DB1405	DB1405	DB24		DB1312/DB1400
CT. 1 DB1405	DB1405	DB24		DB1312/DB1400
CU. 1 DB1405	DB1405	DB24		DB1312/DB1400
CV. 1 DB1405	DB1405	DB24		DB1312/DB1400
CW. 1 DB1405	DB1405	DB24		DB1312/DB1400
CX. 1 DB1405	DB1405	DB24		DB1312/DB1400
CY. 1 DB1405	DB1405	DB24		DB1312/DB1400
CA. 1 DB1405	DB1405	DB24		DB1312/DB1400
CB. 1 DB1405	DB1405	DB24		DB1312/DB1400
CC. 1 DB1405	DB1405	DB24		DB1312/DB1400
CD. 1 DB1405	DB1405	DB24		DB1312/DB1400
CE. 1 DB1405	DB1405	DB24		DB1312/DB1400
CF. 1 DB1405	DB1405	DB24		DB1312/DB1400
CG. 1 DB1405	DB1405	DB24		DB1312/DB1400
CH. 1 DB1405	DB1405	DB24		DB1312/DB1400
CI. 1 DB1405	DB1405	DB24		DB1312/DB1400
CJ. 1 DB1405	DB1405	DB24		DB1312/DB1400
CK. 1 DB1405	DB1405	DB24		DB1312/DB1400
CL. 1 DB1405	DB1405	DB24		DB1312/DB1400
CM. 1 DB1405	DB1405	DB24		DB1312/DB1400
CN. 1 DB1405	DB1405	DB24		DB1312/DB1400
CO. 1 DB1405	DB1405	DB24		DB1312/DB1400
CP. 1 DB1405	DB1405	DB24		DB1312/DB1400
CQ. 1 DB1405	DB1405	DB24		DB1312/DB1400
CR. 1 DB1405	DB1405	DB24		DB1312/DB1400
CS. 1 DB1405	DB1405	DB24		DB1312/DB1400
CT. 1 DB1405	DB1405	DB24		DB1312/DB1400
CU. 1 DB1405	DB1405	DB24		DB1312/DB1400
CV. 1 DB1405	DB1405	DB24		DB1312/DB1400
CW. 1 DB1405	DB1405	DB24		DB1312/DB1400
CX. 1 DB1405	DB1405	DB24		DB1312/DB1400
CY. 1 DB1405	DB1405	DB24		DB1312/DB1400
CA. 1 DB1405	DB1405	DB24		DB1312/DB1400
CB. 1 DB1405	DB1405	DB24		DB1312/DB1400
CC. 1 DB1405	DB1405	DB24		DB1312/DB1400
CD. 1 DB1405	DB1405	DB24		DB1312/DB1400
CE. 1 DB1405	DB1405	DB24		DB1312/DB1400
CF. 1 DB1405	DB1405	DB24		DB1312/DB1400
CG. 1 DB1405	DB1405	DB24		DB1312/DB1400
CH. 1 DB1405	DB1405	DB24		DB1312/DB1400
CI. 1 DB1405	DB1405	DB24		DB1312/DB1400
CJ. 1 DB1405	DB1405	DB24		DB1312/DB1400
CK. 1 DB1405	DB1405	DB24		DB1312/DB1400
CL. 1 DB1405	DB1405	DB24		DB1312/DB1400
CM. 1 DB1405	DB1405	DB24		DB1312/DB1400
CN. 1 DB1405	DB1405	DB24		DB1312/DB1400
CO. 1 DB1405	DB1405	DB24		DB1312/DB1400
CP. 1 DB1405	DB1405	DB24		DB1312/DB1400
CQ. 1 DB1405	DB1405	DB24		DB1312/DB1400
CR. 1 DB1405	DB1405	DB24		DB1312/DB1400
CS. 1 DB1405	DB1405	DB24		DB1312/DB1400
CT. 1 DB1405	DB1405	DB24		DB1312/DB1400
CU. 1 DB1405	DB1405	DB24		DB1312/DB1400
CV. 1 DB1405	DB1405	DB24		DB1312/DB1400
CW. 1 DB1405	DB1405	DB24		DB1312/DB1400
CX. 1 DB1405	DB1405	DB24		DB1312/DB1400
CY. 1 DB1405	DB1405	DB24		DB1312/DB1400
CA. 1 DB1405	DB1405	DB24		DB1312/DB1400
CB. 1 DB1405	DB1405	DB24		DB1312/DB1400
CC. 1 DB1405	DB1405	DB24		DB1312/DB1400
CD. 1 DB1405	DB1405	DB24		DB1312/DB1400
CE. 1 DB1405	DB1405	DB24		DB1312/DB1400
CF. 1 DB1405	DB1405	DB24		DB1312/DB1400
CG. 1 DB1405	DB1405	DB24		DB1312/DB1400
CH. 1 DB1405	DB1405	DB24		DB1312/DB1400
CI. 1 DB1405	DB1405	DB24		DB1312/DB1400
CJ. 1 DB1405	DB1405	DB24		DB1312/DB1400
CK. 1 DB1405	DB1405	DB24		DB1312/DB1400
CL. 1 DB1405	DB1405	DB24		DB1312/DB1400
CM. 1 DB1405	DB1405	DB24		DB1312/DB1400
CN. 1 DB1405	DB1405	DB24		DB1312/DB1400
CO. 1 DB1405	DB1405	DB24		DB1312/DB1400
CP. 1 DB1405	DB1405	DB24		DB1312/DB1400
CQ. 1 DB1405	DB1405	DB24		DB1312/DB1400
CR. 1 DB1405	DB1405	DB24		DB1312/DB1400</

25

(**) Primer DRB22 is used to sequence the negative strand whenever new allelic sequences are identified.

(***) Each DQA1 sequencing primer anneals to a different strand. Reaction 1 uses an alternative amplification primer (DRB122 instead of DRB11) in order to determine the sequence of the negative strand. It may not be expected according to the rest of the haplotype. Reaction 2 is used for sequencing the negative strand of DQB1 in situations where new allelic sequences are identified.

(****) Sequencing of the third exon is necessary to distinguish among certain DPA1 alleles.

(*****) Primer DRB1400 may be used in sequencing amplified DRB genes from genomic DNA to read the sequences immediately following the 3' amplification primer. Reaction 0, 1, 2 and AD are alternative combinations for confirming the sequence of the third exon. It may not be expected according to the rest of the Class II haplotype.

(#) This primer combination is used to distinguish between DRB140701 and DRB140702 which differ by a single base pair in their third exons.

26

EXAMPLE II
Protocol: HLA Class II "Typing" by Direct Sequencing of DRB, DQB, DQA, DPA and DPB Genes

1. Cell Lines and Subjects

5 Lymphoblastoid cell lines (LCLs) representing each of the known Class II haplotypes defined at the 10th International Histocompatibility Workshop [Dupont, Hum. Immunol., 25, 3 (1989)] were provided by Dr. Miriam Segall (University of Minnesota). Forty unrelated subjects who had been previously serologically typed for Class I and Class II antigens were also studied. The serological types of each of the subjects under study were not known to the investigator performing the

sequence analysis at the time the analysis was performed. These subjects included both healthy and affected (insulin-dependent diabetes and autoimmune thyroid disease) individuals. The sequenced haplotypes, many in heterozygote combinations, included: DR7 (n=3),

15 DR17 (n=26), DR4 (n=16), DR11 (n=8), DR8 (n=4), DR1 (n=6), DR15 (n=6), DR16 (n=2), DR13 (n=2), DR14 (n=2), DR212 (n=3), DR526 (n=3). The cell lines and heterozygote combinations tested are shown in Table III.

20 Since the complexity at DPA and DPB loci is similar to that of DQ genes, the primer combinations for DPA and DPB typing were optimized in a smaller group of homozygote and heterozygote subjects.

25

2. HLA-DQB, DQB and DPA Transcript Amplification Using Conserved and Non-Conserved Oligonucleotides

30 Total cellular RNA was prepared from (1 µg) from 5-10x10⁶ peripheral blood mononuclear cells (PBMC) or lymphoblastoid cell lines (LCLs) by cesium chloride centrifugation [Chirgwin et al. Biochemistry, 18, 5249 (1979)]. Alternatively, total RNA from peripheral blood

35 (2-10 ml) was partially purified using a much faster protocol [Gounh, Anal. Biochem., 173, 93 (1988)]. One microgram of total cellular RNA was reverse transcribed with Moloney leukemia virus reverse transcriptase

(MUTVET) (200 u, Bethesda Research Laboratories) in 50 mM Tris HCl, pH 8.3, 75 mM KCl, 10 mM DTT, 3mM MgCl₂, in the presence of the ribonuclease inhibitor RNAs in (10 units, Promega), 75 pH each dNTP and 10 pmole of a specific non-sense primer (Table II) in a 20 ml final volume for 30-45 min at 37°C. Eight μ l of 10X PCR buffer (500 mM KCl, 100 mM Tris-HCl, pH 8.3, 7.5-15 mM MgCl₂, 0.1% gelatin) was added after the incubation period. A 5'-primer (20 pmole) (Type I or Type 2

10 primers, respectively, see Table II) plus 10 pmole of the non-sense primer and two units of Taq polymerase were also added and the final volume was adjusted to 100 μ l with distilled water. The reaction mixture was

15 subjected to 35 cycles of 30 sec at 94°C, 30 sec at 37°C or 55°C and 30 sec at 72°C using a Perkin-Elmer Cetus Thermocycler [see Saiki et al., *Science* (1985); Mullis and Faloona, *Science* (1987); Saiki et al., *Science* (1986);

Scharf et al., *Science* (1986)]. The primers used here,

20 they are shown in Table II. The reactions for each locus are usually performed in separate microfuge tubes. However, when using conserved primers, the cDNA and PCR reactions for all loci (DHB, DQA, DQB, DPA and DPB) can be successfully performed simultaneously in the same tube.

3. Direct Sequencing of Amplified Products with Tag Polymerase

The reaction mixture (100 μ l) was freed of 30 unincorporated dNTPs and excess of oligonucleotides by spin-dialysis using Centricon-100 (Amicon) or Ultrafree-100 (Millipore) microconcentrators. One half of the reaction (20 μ l) was dried down and resuspended in 15 μ l of IX Tag sequencing buffer (50 mM Tris-HCl, pH 9, 10 mM MgCl₂). Internal oligonucleotides were used for 35 priming the sequencing of DQB, DHB, DQA, DPB and DPA genes, respectively (Table II). Primers for sequencing each strand are listed in Table II. Only one strand is

28 routinely sequenced for typing; sequencing of the other strand is performed in cases where a new allelic sequence is suspected. Eighty to 100 ng of primer were end-labelled with 10 pmole of gamma-³²P labelled ATP

5 (5000 Ci/mole, 10 μ Ci/ μ l) and 5 units of T4 polynucleotide kinase (Promega Biotech) in a 10 μ l final volume. Ten ng of primer (1 μ l) were added to the

sequencing mixture without extraction of unincorporated labelled ATP, boiled for 5 min., and then left at room temperature for 15 min. Eight units of recombinant Taq polymerase (USB) were added to the mixture. Four μ l of

the annealed primer/template mixture were later added to 4 μ l of each of the stop nucleotide mixes: a) Term mix

15 ddG: 15 microM each dCTP, dATP, dCTP, dTTP; 45 microM dGTP; b) Term mix dda: 15 microM each dCTP, dATP, dCTP, dTTP; 600 microM ddATP; c) Term mix ddt: 15 microM each dCTP, dATP, dCTP, dTTP; 1200 microM ddCTP; d) Term

mix ddc: 15 microM each dCTP, dATP, dCTP, dTTP; 450 microM ddCTP. The reactions were allowed to proceed for

20 two consecutive periods of 10 min. at 72-74°C. After the second cycle, each reaction was chased with 2 μ l of a 7.5 μ M mixture of ATP, GTP, TTP, CTP, and allowed to proceed for 5 min. After spinning down, the reaction

25 was stopped by adding 4 ml of 95% (vol/vol) formamide/20 mM EDTA, heated to 80°C for 5 min. and loaded on a .4 mm thick 6% polyacrylamide/7M urea gel. Electrophoresis was performed at 2500 V for 2 hr, the gel fixed in 5% (vol/vol) glacial acetic acid/5% (v/v) methanol for 15 min, dried, and exposed to Kodak X-Omat film for 4 to 12

30 hours.

RESULTS

1. Sequence-Based Typing of DR and DQ Polymorphic Genes in Homozygous Typing Cells

35 Homozygous lymphoblastoid cell lines (LCLs)

from the panel of the 10th International Histocompatibility Workshop (Table III) were used as an initial test of the methodology. In total, these cell

lines were representative of most of the known DR and DQ alleles at the time the analysis was conducted.

Total cellular RNA isolated from homozygous

LCMs was reverse-transcribed and the resultant cDNAs amplified using conserved oligonucleotides specific for DRB1/DRB3/DRB4/DRB5 or DQB1 or DQA1 genes as described in the preceding protocol. The conserved or Type 1 oligonucleotide primers anneal to regions of conserved

DNA sequences; these regions are identical among the known alleles at each locus and flank the second exon of

Class II genes. These conserved primers, as opposed to non-conserved or Type 2 primers, are designed to amplify all known alleles at DRB, DQA1 and DQB1 loci and, thus, all possible combinations of these alleles in any given

heterozygote. The Type 1 oligonucleotides did not cross-amplify templates at loci other than those specified by the oligonucleotides (i.e., the DQA1

primers did not amplify DRB or DQB1 transcripts and vice versa); as expected, the DRB primers also amplified any

DRB1, DRB4 or DRB5 transcripts present in addition to DRB1. Sequencing of these amplified templates was

performed using a Type 1 primer annealing to a conserved region of the cDNAs internal to the sequence recognized by the amplification primers. Figure 1A shows the

general strategy for the method (SRT) and Figure 1B shows the relative position of each of the

oligonucleotides used for the cDNA, PCR and sequencing reactions on the mature DRB, DQA and DQB mRNA molecules.

The sequences of these primers, the loci they are

specific for, the specific positions (codons) to which they anneal and the reaction(s) they are used in are

indicated in Table II where the specific combinations of primers that can be used for the cDNA/PCR/sequencing reactions for each locus are identified. As noted in

the legend to Table II, some of the primer combinations shown represent alternatives which may be useful in confirming results for a particular locus which do not

fit with the expected sequences usually found with the rest of the haplotype. Each cDNA/PCR reaction is usually performed in a separate tube. However, when

using Type 1 primers, the cDNA/PCR reactions for all the 5 loci (DRB, DQA, DQB, DPA and DPB) can be performed simultaneously in the same tube. The products of each locus are sequenced in separate tubes. Following the

conditions described in the above protocol, the sequence ladders between the sequencing primer and the 5'

10 amplification primer could be clearly read starting from 2 to 14 bases from the sequencing primer binding site.

No anomalous amplification products or sequencing

ladders were detectable upon direct sequence analysis of amplified DRB, DQB1 and DQA cDNAs from the 43

15 homozygous cell lines tested (Table IIIa). The specific alleles at each Class II HLA locus composing the haplotypes carried by each of these cell lines are shown

in Table IIId. The number of ladders generated for each cell line was always that expected according to the

20 specificity of the amplification primers (one DQB1 and one DQA1 ladder for all cell lines, one DRB ladder for DR1 and DRW8 cell lines and two DRB ladders for

haplotypes of the DRB2 and DRB3 supertypic groups).

Thus, analysis of the homozygous typing cell lines

25 showed that the Type 1 primers used for cDNA synthesis, PCR and

sequencing reactions allowed for accurate amplification and sequencing of all the tested alleles at each of these loci.

TABLE IIIA

Cell Lines and Heterozygote Combinations Tested

Cell Line	Class II HLA Type	Subject	Class II HLA Types**
SA	DR1-Dw1	S1	DR1-Dw1/DRw17
SA70762	DRw15-Dw2	S2	DRw15-Dw2/DRw4
KAS011	DRw16-Dw21	S3	DRw16-Dw21/DRw1
*CALAGHERO	DRw15-Dw1	S4	DRw15-Dw2/DRw17
*K3R076	DRw15-Dw21	S5	DRw15-Dw2/DRw1
*D2X	DRw16-Dw21	S6	DRw16-Dw21/DRw17
WT74	DRw15-Dw21	S7	DRw16-Dw21/DRw17
RXL	DRw15-Dw22	S8	*DRw13-
Dw18/DRw17			
WFO	DRw15-Dw2	S9	DRw13-Dw19/DRw17
WPE	DRw15-Dw2	S10	DRw13-Dw19/DRw17
*AAAI	DRw15-Dw2	S11	DRw13-Dw19/DRw17
E4181324	DRw15-Dw12	S12	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S13	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S14	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S15	*DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S16	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S17	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S18	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S19	*DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S20	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S21	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S22	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S23	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S24	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S25	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S26	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S27	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S28	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S29	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S30	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S31	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S32	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S33	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S34	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S35	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S36	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S37	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S38	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S39	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S40	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S41	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S42	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S43	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S44	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S45	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S46	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S47	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S48	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S49	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S50	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S51	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S52	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S53	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S54	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S55	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S56	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S57	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S58	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S59	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S60	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S61	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S62	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S63	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S64	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S65	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S66	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S67	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S68	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S69	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S70	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S71	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S72	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S73	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S74	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S75	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S76	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S77	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S78	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S79	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S80	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S81	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S82	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S83	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S84	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S85	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S86	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S87	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S88	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S89	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S90	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S91	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S92	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S93	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S94	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S95	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S96	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S97	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S98	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S99	DRw13-Dw19/DRw17
WT148	DRw15-Dw12	S100	DRw13-Dw19/DRw17

The allelic composition at DRB, DQA1 and DQB1 loci for the sequenced haplotypes corresponded to that expected according to published sequence information from well

characterized homozygous cell lines unless indicated (*).

- * Haplotypes carrying new allelic sequences (DRB1, DRB3, DQA1 or DQB1 loci).
- * Only the tested heterozygote combinations are listed. The remainder of the 40 subjects tested were homozygotes or carried the haplotypes listed in this table.
- * This DRB specificity (DRB56) has been given this arbitrary designation according to serological, RFLP and sequence information.

Table VIII

Allelic Composition of Human Class II Haplotypes

Haplotype	DRB1	DRB3	DRB5	DRB1	DRB1
DRB1-Dw1	*0101	-	-	-	-
DRB1-Dw2	*0102	-	-	-	-
DRB1-Dw2.1	*0103	-	-	-	-
DRB1-Dw2.2	*0104	-	-	-	-
DRB1-Dw2.3	*0105	-	-	-	-
DRB1-Dw2.4	*0106	-	-	-	-
DRB1-Dw2.5	*0107	-	-	-	-
DRB1-Dw2.6	*0108	-	-	-	-
DRB1-Dw2.7	*0109	-	-	-	-
DRB1-Dw2.8	*0110	-	-	-	-
DRB1-Dw2.9	*0111	-	-	-	-
DRB1-Dw2.10	*0112	-	-	-	-
DRB1-Dw2.11	*0113	-	-	-	-
DRB1-Dw2.12	*0114	-	-	-	-
DRB1-Dw2.13	*0115	-	-	-	-
DRB1-Dw2.14	*0116	-	-	-	-
DRB1-Dw2.15	*0117	-	-	-	-
DRB1-Dw2.16	*0118	-	-	-	-
DRB1-Dw2.17	*0119	-	-	-	-
DRB1-Dw2.18	*0120	-	-	-	-
DRB1-Dw2.19	*0121	-	-	-	-
DRB1-Dw2.20	*0122	-	-	-	-
DRB1-Dw2.21	*0123	-	-	-	-
DRB1-Dw2.22	*0124	-	-	-	-
DRB1-Dw2.23	*0125	-	-	-	-
DRB1-Dw2.24	*0126	-	-	-	-
DRB1-Dw2.25	*0127	-	-	-	-
DRB1-Dw2.26	*0128	-	-	-	-
DRB1-Dw2.27	*0129	-	-	-	-
DRB1-Dw2.28	*0130	-	-	-	-
DRB1-Dw2.29	*0131	-	-	-	-
DRB1-Dw2.30	*0132	-	-	-	-
DRB1-Dw2.31	*0133	-	-	-	-
DRB1-Dw2.32	*0134	-	-	-	-
DRB1-Dw2.33	*0135	-	-	-	-
DRB1-Dw2.34	*0136	-	-	-	-
DRB1-Dw2.35	*0137	-	-	-	-
DRB1-Dw2.36	*0138	-	-	-	-
DRB1-Dw2.37	*0139	-	-	-	-
DRB1-Dw2.38	*0140	-	-	-	-
DRB1-Dw2.39	*0141	-	-	-	-
DRB1-Dw2.40	*0142	-	-	-	-
DRB1-Dw2.41	*0143	-	-	-	-
DRB1-Dw2.42	*0144	-	-	-	-
DRB1-Dw2.43	*0145	-	-	-	-
DRB1-Dw2.44	*0146	-	-	-	-
DRB1-Dw2.45	*0147	-	-	-	-
DRB1-Dw2.46	*0148	-	-	-	-
DRB1-Dw2.47	*0149	-	-	-	-
DRB1-Dw2.48	*0150	-	-	-	-
DRB1-Dw2.49	*0151	-	-	-	-
DRB1-Dw2.50	*0152	-	-	-	-
DRB1-Dw2.51	*0153	-	-	-	-
DRB1-Dw2.52	*0154	-	-	-	-
DRB1-Dw2.53	*0155	-	-	-	-
DRB1-Dw2.54	*0156	-	-	-	-
DRB1-Dw2.55	*0157	-	-	-	-
DRB1-Dw2.56	*0158	-	-	-	-
DRB1-Dw2.57	*0159	-	-	-	-
DRB1-Dw2.58	*0160	-	-	-	-
DRB1-Dw2.59	*0161	-	-	-	-
DRB1-Dw2.60	*0162	-	-	-	-
DRB1-Dw2.61	*0163	-	-	-	-
DRB1-Dw2.62	*0164	-	-	-	-
DRB1-Dw2.63	*0165	-	-	-	-
DRB1-Dw2.64	*0166	-	-	-	-
DRB1-Dw2.65	*0167	-	-	-	-
DRB1-Dw2.66	*0168	-	-	-	-
DRB1-Dw2.67	*0169	-	-	-	-
DRB1-Dw2.68	*0170	-	-	-	-
DRB1-Dw2.69	*0171	-	-	-	-
DRB1-Dw2.70	*0172	-	-	-	-
DRB1-Dw2.71	*0173	-	-	-	-
DRB1-Dw2.72	*0174	-	-	-	-
DRB1-Dw2.73	*0175	-	-	-	-
DRB1-Dw2.74	*0176	-	-	-	-
DRB1-Dw2.75	*0177	-	-	-	-
DRB1-Dw2.76	*0178	-	-	-	-
DRB1-Dw2.77	*0179	-	-	-	-
DRB1-Dw2.78	*0180	-	-	-	-
DRB1-Dw2.79	*0181	-	-	-	-
DRB1-Dw2.80	*0182	-	-	-	-
DRB1-Dw2.81	*0183	-	-	-	-
DRB1-Dw2.82	*0184	-	-	-	-
DRB1-Dw2.83	*0185	-	-	-	-
DRB1-Dw2.84	*0186	-	-	-	-
DRB1-Dw2.85	*0187	-	-	-	-
DRB1-Dw2.86	*0188	-	-	-	-
DRB1-Dw2.87	*0189	-	-	-	-
DRB1-Dw2.88	*0190	-	-	-	-
DRB1-Dw2.89	*0191	-	-	-	-
DRB1-Dw2.90	*0192	-	-	-	-
DRB1-Dw2.91	*0193	-	-	-	-
DRB1-Dw2.92	*0194	-	-	-	-
DRB1-Dw2.93	*0195	-	-	-	-
DRB1-Dw2.94	*0196	-	-	-	-
DRB1-Dw2.95	*0197	-	-	-	-
DRB1-Dw2.96	*0198	-	-	-	-
DRB1-Dw2.97	*0199	-	-	-	-
DRB1-Dw2.98	*0200	-	-	-	-
DRB1-Dw2.99	*0201	-	-	-	-
DRB1-Dw2.100	*0202	-	-	-	-
DRB1-Dw2.101	*0203	-	-	-	-
DRB1-Dw2.102	*0204	-	-	-	-
DRB1-Dw2.103	*0205	-	-	-	-
DRB1-Dw2.104	*0206	-	-	-	-
DRB1-Dw2.105	*0207	-	-	-	-
DRB1-Dw2.106	*0208	-	-	-	-
DRB1-Dw2.107	*0209	-	-	-	-
DRB1-Dw2.108	*0210	-	-	-	-
DRB1-Dw2.109	*0211	-	-	-	-
DRB1-Dw2.110	*0212	-	-	-	-
DRB1-Dw2.111	*0213	-	-	-	-
DRB1-Dw2.112	*0214	-	-	-	-
DRB1-Dw2.113	*0215	-	-	-	-
DRB1-Dw2.114	*0216	-	-	-	-
DRB1-Dw2.115	*0217	-	-	-	-
DRB1-Dw2.116	*0218	-	-	-	-
DRB1-Dw2.117	*0219	-	-	-	-
DRB1-Dw2.118	*0220	-	-	-	-
DRB1-Dw2.119	*0221	-	-	-	-
DRB1-Dw2.120	*0222	-	-	-	-
DRB1-Dw2.121	*0223	-	-	-	-
DRB1-Dw2.122	*0224	-	-	-	-
DRB1-Dw2.123	*0225	-	-	-	-
DRB1-Dw2.124	*0226	-	-	-	-
DRB1-Dw2.125	*0227	-	-	-	-
DRB1-Dw2.126	*0228	-	-	-	-
DRB1-Dw2.127	*0229	-	-	-	-
DRB1-Dw2.128	*0230	-	-	-	-
DRB1-Dw2.129	*0231	-	-	-	-
DRB1-Dw2.130	*0232	-	-	-	-
DRB1-Dw2.131	*0233	-	-	-	-
DRB1-Dw2.132	*0234	-	-	-	-
DRB1-Dw2.133	*0235	-	-	-	-
DRB1-Dw2.134	*0236	-	-	-	-
DRB1-Dw2.135	*0237	-	-	-	-
DRB1-Dw2.136	*0238	-	-	-	-
DRB1-Dw2.137	*0239	-	-	-	-
DRB1-Dw2.138	*0240	-	-	-	-
DRB1-Dw2.139	*0241	-	-	-	-
DRB1-Dw2.140	*0242	-	-	-	-
DRB1-Dw2.141	*0243	-	-	-	-
DRB1-Dw2.142	*0244	-	-	-	-
DRB1-Dw2.143	*0245	-	-	-	-
DRB1-Dw2.144	*0246	-	-	-	-
DRB1-Dw2.145	*0247	-	-	-	-
DRB1-Dw2.146	*0248	-	-	-	-
DRB1-Dw2.147	*0249	-	-	-	-
DRB1-Dw2.148	*0250	-	-	-	-
DRB1-Dw2.149	*0251	-	-	-	-
DRB1-Dw2.150	*0252	-	-	-	-
DRB1-Dw2.151	*0253	-	-	-	-
DRB1-Dw2.152	*0254	-	-	-	-
DRB1-Dw2.153	*0255	-	-	-	-
DRB1-Dw2.154	*0256	-	-	-	-
DRB1-Dw2.155	*0257	-	-	-	-
DRB1-Dw2.156	*0258	-	-	-	-
DRB1-Dw2.157	*0259	-	-	-	-
DRB1-Dw2.158	*0260	-	-	-	-
DRB1-Dw2.159	*0261	-	-	-	-
DRB1-Dw2.160	*0262	-	-	-	-
DRB1-Dw2.161	*0263	-	-	-	-
DRB1-Dw2.162	*0264	-	-	-	-
DRB1-Dw2.163	*0265	-	-	-	-
DRB1-Dw2.164	*0266	-	-	-	-
DRB1-Dw2.165	*0267	-	-	-	-
DRB1-Dw2.166	*0268	-	-	-	-
DRB1-Dw2.167	*0269	-	-	-	-
DRB1-Dw2.168	*0270	-	-	-	-
DRB1-Dw2.169	*0271	-	-	-	-
DRB1-Dw2.170	*0272	-	-	-	-
DRB1-Dw2.171	*0273	-	-	-	-
DRB1-Dw2.172	*0274	-	-	-	-
DRB1-Dw2.173	*0275	-	-	-	-
DRB1-Dw2.174	*0276	-	-	-	-
DRB1-Dw2.175	*0277	-	-	-	-
DRB1-Dw2.176	*0278	-	-	-	-
DRB1-Dw2.177	*0279	-	-	-	-
DRB1-Dw2.178	*0280	-	-	-	-
DRB1-Dw2.179	*0281	-	-	-	-
DRB1-Dw2.180	*0282	-	-	-	-
DRB1-Dw2.181	*0283	-	-	-	-
DRB1-Dw2.182	*0284	-	-	-	-
DRB1-Dw2.183	*0285	-	-	-	-
DRB1-Dw2.184	*0286	-	-	-	-
DRB1-Dw2.185	*0287	-	-	-	-
DRB1-Dw2.186	*0288	-	-	-	-
DRB1-Dw2.187	*0289	-	-	-	-
DRB1-Dw2.188	*0290	-	-	-	-
DRB1-Dw2.189	*0291	-	-	-	-
DRB1-Dw2.190	*0292	-	-	-	-
DRB1-Dw2.191	*0293	-	-	-	-
DRB1-Dw2.192	*0294	-	-	-	-
DRB1-Dw2.193	*0295	-	-	-	-
DRB1-Dw2.194	*0296	-	-	-	-
DRB1-Dw2.195	*0297	-	-	-	-
DRB1-Dw2.196	*0298	-	-	-	-
DRB1-Dw2.197	*0299	-	-	-	-
DRB1-Dw2.198	*0300	-	-	-	-
DRB1-Dw2.199	*0301	-	-	-	-
DRB1-Dw2.200	*0302	-	-	-	-
DRB1-Dw2.201	*0303	-	-	-	-
DRB1-Dw2.202	*0304	-	-	-	-
DRB1-Dw2.203	*0305	-	-	-	-
DRB1-Dw2.204	*0306	-	-	-	-
DRB1-Dw2.205	*0307	-	-	-	-
DRB1-Dw2.206	*0308	-	-	-	-
DRB1-Dw2.207	*0309	-	-	-	-
DRB1-Dw2.208	*0310	-	-	-	-
DRB1-Dw2.209	*0311	-	-	-	-
DRB1-Dw2.210	*0312	-	-	-	-
DRB1-Dw2.211	*0313	-	-	-	-
DRB1-Dw2.212	*0314	-	-	-	-
DRB1-Dw2.213	*0315	-	-	-	-
DRB1-Dw2.214	*0316	-	-	-	-
DRB1-Dw2.215	*0317	-	-	-	-
DRB1-Dw2.216	*0318	-	-	-	-
DRB1-Dw2.217	*0319	-	-	-	-
DRB1-Dw2.218	*0320	-	-	-	-
DRB1-Dw2.219	*0321	-	-	-	-
DRB1-Dw2.220	*0322	-	-	-	-
DRB1-Dw2.221	*0323	-	-	-	-
DRB1-Dw2.222	*0324	-	-	-	-
DRB1-Dw2.223	*0325	-	-	-	-
DRB1-Dw2.224	*0326	-	-	-	-
DRB1-Dw2.225	*0327	-	-	-	-
DRB1-Dw2.226	*0328	-	-	-	-
DRB1-Dw2.227	*0329	-	-	-	-
DRB1-Dw2.228	*0330	-	-	-	-
DRB1-Dw2.229	*0331	-	-	-	-
DRB1-Dw2.230	*0332	-	-	-	-
DRB1-Dw2.231	*0333	-	-	-	-
DRB1-Dw2.232	*0334	-	-	-	-
DRB1-Dw2.233	*0335	-	-	-	-
DRB1-Dw2.234	*0336	-	-	-	-
DRB1-Dw2.235	*0337	-	-	-	-
DRB1-Dw2.236	*0338	-	-	-	-
DRB1-Dw2.237	*0339	-	-	-	-
DRB1-Dw2.238	*0340	-	-	-	-
DRB1-Dw2.239	*0341	-	-	-	-
DRB1-Dw2.240	*0342	-	-	-	-
DRB1-Dw2.241	*0343	-	-	-	-

2. Amplification and Direct Sequencing of DQAI and DQB1 cDNAs in Subjects of Unknown HLA Type

DNA sequences have been determined for most HLA Class II allelic specificities defined by conventional HLA typing techniques (March, S.G.B., Bodmer, J.G. HLA-DNA nucleotide sequences, 1990. Immunogenetics 31:141, 1990; Todd, J.A., Bell, J.I., McDevitt, H.O.: HLA-DQB genes contribute to susceptibility and resistance to insulin-dependent diabetes mellitus. Nature 329:599, 1987). Comparisons of these sequences indicates that any given DQAI or DQB1 homozygous or heterozygous allelic combination is characterized by a specific sequencing ladder.

Total RNA from PBMCs from 40 different subjects was tested to evaluate if the allelic composition of DQAI and DQB1 homo- and heterozygotes could be determined correctly by direct amplification and sequencing using type 1 primers. These subjects had been previously serologically typed but the typing information was not known to the investigator who assigned the Class II allelic specificities from the sequencing results. These 40 subjects comprised 27 different heterozygote combinations (Table III). All individuals were assigned DQAI and DQB1 allelic sequences that were consistent with the serological phenotypes. In all the heterozygotes tested, both allelic sequences could be read clearly from the composite sequence pattern. A unique pattern is found for every particular heterozygote combination in the same way that certain RFLP banding patterns correspond to certain heterozygote allelic combinations. For instance, in a DQB2/DQB1.1 heterozygote one would find the sequence GGAG(A/T)T(G/C)CGGC(A/G) at codons 45 to 49 which can only be attributed to that particular allele combination. In practice, interpretation of heterozygous sequence ladders is initiated by reading certain polymorphic positions where allele-specific bases may be found, such as, for instance, the second

base of codon 46, where DQB1*0201 is the only allele that has an A. The sequences of the two possible templates are then deduced and compared with the sequences of all known alleles at the different loci.

In Figure 3 we show the overlapping ladder corresponding to a DQB1*0201/DQB1*0302 heterozygote; interpretation of the pattern is indicated on the side of the ladder.

The absence of expected bands or the presence of unexpected bands for a particular allele or allelic combination is therefore suggestive of sequence heterogeneity, i.e., new alleles. The same can be said for DQAI and DQB1 typing when appropriate primer combinations are used (Table III). For instance, substitution of the A at the second base of codon 46 would strongly suggest the presence of a sequence variant of DQB1*0201. Once detected, the sequence of the variant can be confirmed after selective amplification of the variant or by subcloning the amplified products.

3. Amplification and Direct Sequencing of DRB cDNAs from Subjects of Unknown HLA Type

As described above, the use of type 1 primers allows the unambiguous sequencing of all heterozygous combinations of DQAI and DQB1 alleles. The same can be said for DPAI and DPBI typing when appropriate primer combinations are used (Table II). Because of the isotypic complexity of DRB genes (expression of more than one DRB locus by certain haplotypes), amplification and sequencing of cDNAs from DRB heterozygotes with type 1 primers can generate up to four overlapping ladders, thus generating complex sequencing patterns.

DRB cDNAs from the same 40 individual tested above for DQAI and DQB1 genes were amplified and sequenced using DRB-specific type 1 primers. As mentioned above, these 40 individuals comprised 27 different heterozygote combinations, including several examples from each of the groups of complex DRB allelic

37

combinations which would generate up to four sequencing ladders. The DRB sequence ladders generated with Type 1 primers were analyzed as described above for DQA1 and DQB1 loci: highly polymorphic positions were analyzed first for the presence of bands unique to specific alleles or groups of alleles (i.e., DR4) and the sequences deduced and compared with the sequences of all known alleles at all loci. As example, in Figure 3 we show the ladder generated by sequencing a complex DRB heterozygote (four overlapping ladders); the positions with two or more bands are indicated on the side of the Figure and assigned to each of the allelic types composing the complex sequencing pattern. For all but one sample, the information deduced from these sequencing experiments matched the independently determined serological phenotypes of the subject under study as well as the DQA1 and DQB1 allelic types assigned to these individuals by direct sequencing of these genes as described above. The inconsistent sample 20 had been serologically typed as DRW13/DR4 but was typed by sequence analysis as DRW13/DRB8-Dw8.1. The presence of a DRB1-0801 allele instead of a DRB1-0401 allele was confirmed in a repeated experiment; we thus believe that the serological typing was in error. In all the 40 cases, all DQB1, DQA1 and DRB1 templates had been equally amplified and sequenced with a similar efficiency by the use of Type 1 primers. DRB3, DRB4 and DRB5 sequence ladders could be read in all but one case (a DRB3*0101 [DRB*52a] sequence was not initially observed in a DRW13/DRB17 heterozygote). Since DRB3*0101 is in linkage disequilibrium with DRB1*0301, the former allele was expected to be found in the overlapping ladder as well. In order to rule out the possibility of an error, the investigator assigning the 35 HLA types from the sequencing ladders repeated the typing of this individual; the DRB3*0101 could be read in the repeated experiment.

WO 92/15711

PCT/US92/01675

38

Although the results generated by the use of Type 1 primers were compatible with the serological phenotypes, the exclusive use of Type 1 primers will not allow in all cases to assign each of the specific 5 ladders to each of the expressed loci in all possible heterozygotes. Given below are the most complex situations which cannot be addressed by the exclusive use of Type 1 primers: 1) distinction among the different DR4 allelic sequences in certain heterozygotes since they differ by only a few nucleotide base pairs and such differences could be masked by the presence of additional ladders; 2) to distinguish between DRB1*1601 and DRB1*1502 since their sequence differences will be masked by those of their linked DRB5 alleles; 3) to distinguish between DRB1*1301 and DRB1*1302 (which only differ at codon 86 since this difference can also be masked by other ladders; and finally 4) distinction between DRB1*0301 and DRB1*0302 in specific heterozygote combinations.

We have thus developed a more informative strategy to deal with DRB; this strategy, which consists of the additional use of non-conserved (Type 2) primers permits the clear elucidation of even the most complex combination of the four DRB sequences that might be present in an individual. These non-conserved primers, as opposed to allele-specific primers, are designed to be used in reactions performed simultaneously with the reactions using Type 1 primers and aim at selectively amplifying certain ladders from the complex sequencing patterns without requiring previous typing information.

Analysis of the sequence variability of the second exon of the DRB genes has allowed us to identify two regions which could be used to design non-conserved (Type 2) primers: 1) codons 5-13; and 2) codons 29-35. The sequence of the former region follows a group-specific sequence pattern, i.e., a sequence shared by groups of alleles at individual loci. The latter region

39

exhibits a scattered nucleotide polymorphism in DRB1 and DRB3, DRB4 and DRB5 genes. We designed five different non-conserved primers annealing to these two polymorphic regions: 1) DRB23 (specific for DR2-DRB1 ladders); 2)

5 DRB24 (specific for DRW17-, DRW18-, DRW13-, DRW14-, DRW11-, DRW12-, and DRW8- DRB1 ladders); 3) DRB25

(specific for DR4- DRB1 ladders); 4) DRB16 and 5) DRB17, the latter two primers annealing to the second region of

moderate polymorphism (from 1 to 5 nucleotides different among the known alleles for each locus) (Tables I and

10 IV). Because of the different nature and distribution of mismatches between these primers and the different DRB templates, the type of templates amplified

selectively by these primers will be different. Each of the first three primers will amplify up to two DRB1

15 cDNAs in any given heterozygote and will not amplify any DRB3, DRB4 or DRB5 cDNAs. On the contrary, the use of primers DRB16 and DRB17 will allow the random selective amplification of certain transcripts from DRB1, DRB3,

20 DRB4 and/or DRB5 loci in most heterozygote combinations. We therefore tested these primers in order to determine which combination would give the best discriminatory

25 results for DRB typing. Furthermore, since the sequences of these primers carry from 0-12 mismatches with the sequences of the known DRB alleles at the

different DRB loci, their use allowed us to determine the number of mismatches between the primers and each of the possible cDNAs that are required to obtain such

selective amplification of DRB transcripts. The specific combinations of primers used for the

30 cDNA/PCR/sequencing reactions are shown in Table II above. The results of this analysis are shown below.

Table II. Mismatches between Type 2 DRB primers and DRB alleles at different loci.

Contribution of Nucleotide Base Pair Mismatches between 5' Amplification Primers and DRB Alleles to the Selective Amplification of Allelic

and/or Non-Allelic DRB Transcripts

TABLE IV

DRB1	DRB3/DRB4/DRB5	DRB3*0101	DRB3*0201	DRB3*0301	DRB4
DRB16	0	4	5	4	2
DRB17	4	3	1	5	2
DRB24	0	6	4	5	4
DRB25	4	3	4	5	4
DRB26	4	3	4	5	4
DRB27	4	3	4	5	4
DRB28	4	3	4	5	4
DRB29	4	3	4	5	4
DRB30	4	3	4	5	4
DRB31	4	3	4	5	4
DRB32	4	3	4	5	4
DRB33	4	3	4	5	4
DRB34	4	3	4	5	4
DRB35	4	3	4	5	4
DRB36	4	3	4	5	4
DRB37	4	3	4	5	4
DRB38	4	3	4	5	4
DRB39	4	3	4	5	4
DRB40	4	3	4	5	4
DRB41	4	3	4	5	4
DRB42	4	3	4	5	4
DRB43	4	3	4	5	4
DRB44	4	3	4	5	4
DRB45	4	3	4	5	4
DRB46	4	3	4	5	4
DRB47	4	3	4	5	4
DRB48	4	3	4	5	4
DRB49	4	3	4	5	4
DRB50	4	3	4	5	4
DRB51	4	3	4	5	4
DRB52	4	3	4	5	4
DRB53	4	3	4	5	4
DRB54	4	3	4	5	4
DRB55	4	3	4	5	4
DRB56	4	3	4	5	4
DRB57	4	3	4	5	4
DRB58	4	3	4	5	4
DRB59	4	3	4	5	4
DRB60	4	3	4	5	4
DRB61	4	3	4	5	4
DRB62	4	3	4	5	4
DRB63	4	3	4	5	4
DRB64	4	3	4	5	4
DRB65	4	3	4	5	4
DRB66	4	3	4	5	4
DRB67	4	3	4	5	4
DRB68	4	3	4	5	4
DRB69	4	3	4	5	4
DRB70	4	3	4	5	4
DRB71	4	3	4	5	4
DRB72	4	3	4	5	4
DRB73	4	3	4	5	4
DRB74	4	3	4	5	4
DRB75	4	3	4	5	4
DRB76	4	3	4	5	4
DRB77	4	3	4	5	4
DRB78	4	3	4	5	4
DRB79	4	3	4	5	4
DRB80	4	3	4	5	4
DRB81	4	3	4	5	4
DRB82	4	3	4	5	4
DRB83	4	3	4	5	4
DRB84	4	3	4	5	4
DRB85	4	3	4	5	4
DRB86	4	3	4	5	4
DRB87	4	3	4	5	4
DRB88	4	3	4	5	4
DRB89	4	3	4	5	4
DRB90	4	3	4	5	4
DRB91	4	3	4	5	4
DRB92	4	3	4	5	4
DRB93	4	3	4	5	4
DRB94	4	3	4	5	4
DRB95	4	3	4	5	4
DRB96	4	3	4	5	4
DRB97	4	3	4	5	4
DRB98	4	3	4	5	4
DRB99	4	3	4	5	4
DRB100	4	3	4	5	4

DRB5 gene from cell line AHA1 has an additional nucleotide substitution in the first base of codon 30, in comparison with DRB5 genes of other DR2 haplotypes. # The DRB1 gene of this specificity (DRB56) has been given this arbitrary designation according to serological, RFLP and sequence information.

The ability of non-conserved primers to select certain alleles in heterozygote combinations was also tested for DQB1 genes (Table VI). As with DRB-specific type 2 primers, the use of high temperatures (55°C) in the annealing step of the PCR was required for achieving the selective amplification of single DQB1 alleles in heterozygotes with non-conserved primers. For instance, when annealing of primer DQB5 was allowed to proceed at 37°C in cDNAs from a DQB1*0301/DQB1*0501 and a DQB1*0201/DQB1*0603 heterozygote, both alleles in both heterozygotes were equally amplified. At 55°C, the allele with the most homologous sequence to the 5' primer, was amplified over the other in the PCR. Combinations of alleles both differing from the primer in two nucleotides but in different relative position were also differentially amplified with a non-conserved primer. For instance, primer DQB5 selected the DQB1*0604 sequence in a DQB1*0604/DQB1*0502 heterozygote (Table VI). Five nucleotides separate the two mismatches between the DQB1*0604 allele and the DQB5 primer, whereas only two nucleotides separate the 20 mismatches between the DQB1*0502 and the primer.

These results clearly indicate that the oligonucleotide primers annealing to polymorphic regions at the 5' end of the target cDNAs can be tailored to achieve a reproducible selective amplification of a limited number of DRB or DQB templates in complex heterozygous combinations. Although the use of Type 1 primers allows the unambiguous sequencing of all possible DQA, DQB, DPA and DRB heterozygotes, such an approach will not give absolute discriminatory information for all DRB heterozygotes. We have shown that the simultaneous use of Type 1 and Type 2 primers for DRB will permit the clear elucidation of even the most complex of all DRB heterozygote combinations. When DRB-Sr is used for typing purposes, we perform three Type 2-reactions (using DRB23, 24 and 25) simultaneously with a Type 1-reaction (Table II). The simultaneous use of these

reactions using these primers has the highest discriminatory power for complete DRB typing in a single run and allows the identification of novel sequence heterogeneity. Only one Type 1 reaction is required for DQB1, one for DQA1, one for DPA1 and one for DRB1 (Table II).

EXAMPLE III

Determining Unknown HLA Types of Subjects by Direct Sequencing of the Second Exon of Class II Genes

Routine HLA typing of large populations of individuals for sequence polymorphisms can be performed by the use of the methodology reported here which can also identify previously unknown allelic variants. Figures 2A and 2B show a flow-chart for the routine protocol used to determine sequence allelism of individuals of unknown HLA types.

1. Employment of Primer Combinations for cDNA, PCR and Direct Sequencing Using RNA as Initial Template

For synthesizing cDNA molecules, the present invention provides single strand DNA anti-sense oligonucleotide primers that anneal to conserved regions of the gene mRNAs to be reverse transcribed, amplified and sequenced. These oligonucleotide primers include an oligonucleotide sequence that: (1) anneals to a conserved region (codons 105 through 111) shared by all the alleles at all the DRB loci, the latter being DRB1, DRB3, DRB4 and DRB5, respectively (e.g., primer DRB20). Four simultaneous cDNA reactions (one per tube) are performed for DRB typing, all using primer DRB20 (reactions A, B, C and D in Table II and Figure 2A); (2) anneals to a conserved region (codons 105 through 111) shared by all the alleles at the DQB locus (e.g., primer DQB27) (reaction E in Table II and Figure 2A); (3) anneals to a conserved region (codons 147 through 157) shared by all the alleles at the DPA locus (e.g., primer DPA9)

(reaction F in Table II and Figure 2A); (4) anneals to a conserved region (codons 105 through 111) shared by all the alleles at the DBP locus (e.g. primer DBP11) (reaction G in Table II and Figure 2A); (5) anneals to a conserved region (codons 104 through 110) shared by all the alleles at the DBP locus (e.g. primer DBP14) (reaction H in Table II and Figure 2A); (6) anneals to a conserved region (codons 222 through 228) shared by all the alleles at the DBP locus (e.g. primer DBP19) (reaction I in Table II and Figure 2A). The specific oligonucleotides added to each of these reactions once the cDNA synthesis is done in order to amplify and sequence the products are indicated below as well as in Table II and in Figure 2.

To amplify cDNA molecules corresponding to each expressed DBP loci of each chromosome (DBP1 and DBP3 or DBP4 or DBP5, depending on the haplotype -isotypic complexity-), a conserved oligonucleotide primer which anneals to codons -32 to -26 (e.g. oligonucleotide DBP11) is added to one of the four tubes where the cDNA synthesis reactions corresponding to DBP genes took place. The combination of the cDNA synthesis reaction primer and the newly added conserved primer is used to amplify all the alleles at all DBP loci expressed by a given individual. Each of the remaining three tubes containing DBP cDNA products receives one of three different non-conserved oligonucleotides (also called Type 2) annealing to codons 7-13 (e.g. primer DBP23), 5-11 (e.g. primer DBP24), 6-13 (e.g. primer DBP25), respectively. Each non-conserved primer is designed to favor the amplification of cDNAs corresponding to different groups of alleles at the DBP1 locus. Comparison of the sequencing ladders generated by these four reactions allows complete and accurate interpretation of the sequences corresponding to each of the four possible DBP genes expressed by a given individual (one or two for each of the parental chromosomes).

For the DBP1 locus, a conserved oligonucleotide primer which anneals to codons 1-7 of the DBP cDNAs (e.g. primer DBP13) can be used for amplifying each of the DBP1 genes expressed in any given individual (one for each parental chromosome). In the case of the DBP1 locus, a conserved single strand DNA oligonucleotide primer useful for amplifying each of the DBP1 genes expressed in any given subject anneals to codons -10 to -4 of the DBP1 cDNA (e.g. primer DBP10). For the DBP1 locus, a conserved oligonucleotide (e.g. primer DBP10) annealing to codons -19 to -13, is used to amplify each of the expressed DBP1 genes in any given subject. For DBP1 locus, a conserved oligonucleotide (e.g. primer DBP15) annealing to codons -23 to -17, is used to amplify each of the expressed DBP1 genes in a given subject. In a separate reaction, conserved primer DBP18, annealing to codons 59-65 of the DBP1 cDNAs is used in combination with the cDNA primer DBP19 to amplify each of the expressed DBP1 genes in any individual. This second DBP1 reaction is targeted at a second polymorphic region of this gene.

Primers useful in direct sequencing the polymerase chain reaction products corresponding to DBP loci include an anti-sense oligonucleotide primer (e.g. DBP12) annealing to codons 87-94 of all the alleles at DBP loci; this primer is used for sequencing the products generated by the first of the four DBP reactions. For direct sequencing the polymerase chain reaction products generated with the other three DBP reactions, an anti-sense oligonucleotide annealing to codons 97-103 of all the alleles at DBP1 locus can be used (e.g. primer DBP30). The use of a different sequencing oligonucleotide in these three DBP reactions allows reading of downstream polymorphic regions of DBP1 genes not seen in the first DBP reaction which uses the example sequencing primer DBP12. Primers useful in direct sequencing the polymerase chain reaction products corresponding to DBP1

47

locus include an anti-sense oligonucleotide primer (e.g. DQ85) annealing to codons 78-83 of all the alleles at this locus. Direct sequencing of polymerase chain reaction products corresponding to DQ1 locus include an anti-sense oligonucleotide primer (e.g. DQ429) annealing to codons 88-95 of all the alleles at this locus. Direct sequencing of polymerase chain reaction products corresponding to DPA1 locus include a sense oligonucleotide primer (e.g. DPA13) annealing to codons 12/-5 of all the alleles at this locus.

10 For direct sequencing of polymerase chain reaction products for the DPA1 reaction which used primers DPA14 and DPA15, an anti-sense oligonucleotide annealing to codons 88-94 of all the alleles at this locus can be used (e.g. primer DPA16).

15 For direct sequencing of polymerase chain reaction products for the DPA1 reaction which used primers DPA19 and DPA16, an anti-sense oligonucleotide annealing to codons 214-220 of all the alleles at this locus can be used (e.g. primer DPA20).

2. Employment of Primer Combinations for PCR and Direct Sequencing Using DNA Templates

To amplify DNA molecules corresponding to each DRB loci of each chromosome a conserved anti-sense oligonucleotide primer annealing to base pairs 18-38 of intron 3 (e.g. oligonucleotide DRB1406) is added to each of 25 four PCR reaction tubes (reactions S, V, T and U in Table II and Figure 26). Each of these four tubes will receive a different additional oligonucleotide annealing to codons -4 to +3 (e.g. primer DRB22), to codons 7-13 (e.g. primer DRB23), 5-11 (e.g. primer DRB24), 6-13 (e.g. primer DRB25), respectively. The first reaction is used to amplify all the alleles at all DRB loci carried by a given individual. Each of the remaining three reactions is designed to favor the amplification of DNA corresponding to different groups of alleles at the DRB1 locus. As with RNA templates, comparison 35 of the sequencing ladders generated by these four reactions

48

allows complete and accurate interpretation of the sequences corresponding to each of the four possible DRB genes expressed by a given individual (one or two for each of the parental chromosomes).

5 For the DQB1 locus, two conserved oligonucleotide primers which anneal to codons 88-94 (e.g. primer DQB932) and 11-17 (e.g. primer DQB931) or 1-7 (e.g. primer DQB13) can be used for amplifying each of the DQB1 genes carried by any given individual (one for each parental chromosome) (reaction W in Table II and Figure 2B). For the DPA1 locus, two conserved oligonucleotides (a primer, e.g. DPA14, annealing to base pairs -42 to -62 of intron 2, and a primer e.g. DPA15, annealing base pairs 39-49 of intron 3) are used to amplify each of the DPA1 genes carried by any given subject (reaction X in Table II and Figure 2B). For DPA1 locus a conserved oligonucleotide such as DPA10 (annealing to base pairs -69 to -50 of intron 2) and DPA11 (annealing to base pairs 55-71 of intron 3) are used to amplify each of the DPA1 genes carried by a given subject (reaction Y in Table II and Figure 2B).

20 Primers useful in direct sequencing the polymerase chain reaction products generated from DNA templates corresponding to DRB loci include an anti-sense oligonucleotide primer (e.g. DRB12) annealing to codons 87-94 of all alleles at DRB loci; this primer is used for sequencing the products generated by the first of the four DRB reactions. For direct sequencing the polymerase chain reaction products generated with the other three DRB reactions, a sense oligonucleotide annealing to codons 39-46 of all the alleles at DRB1 locus can be used (e.g. primer DRB1400). The use of a different sequencing oligonucleotide in these three DRB reactions allows reading of downstream polymorphic regions of DRB1 genes not seen in the first DRB reaction which uses the example sequencing primer DRB12.

35 Primers useful in direct sequencing the polymerase chain

reaction products corresponding to DQB1 locus include an anti-sense oligonucleotide primer (e.g. DQB5) annealing to codons 78-83 of all the alleles at this locus. Direct sequencing of polymerase chain reaction products corresponding to DPA1 locus include an anti-sense oligonucleotide primer (e.g. DPA16) annealing to base pairs 1-21 of intron 3 of all the alleles at this locus. For direct sequencing of polymerase chain reaction products for the DPA1 reaction an anti-sense oligonucleotide annealing to codons 76-82 of all the alleles at this locus can be used (e.g. primer DPA12).

Procedure for Determining Unknown HLA Type

A subject of unknown HLA type, diseased or not, is to be typed for Class II HLA polymorphisms. From 10 to 50 μ l of peripheral blood are drawn. The peripheral blood mononuclear cells are prepared by centrifugation over Ficoll-Hypaque gradients. The cells are then lysed in guanidium isothiocyanate and total cellular RNA prepared using conventional methods (either by centrifugation on cesium chloride gradients, which lasts about 16 hours, or by the guanidium isothiocyanate-phenol-chloroform extraction method, which can be performed in less than 4 hours. See Gomb, *et al.* (1988); Johns *et al.*, *Anal. Biochem.* 180:276 (1989). Otherwise genomic DNA from these cells or other sources (hair, blood stains, sperm, etc.) can be prepared with conventional methods such as provided by Higuchi, R. in *PCR Technology*, Erlich, M. (ed.), Stockton Press:131 (1989). DQB1, DQA1, DQB, DQB3, DQB3/4/5, DPA1 and DPA1 cDNA molecules are synthesized from total RNA using locus-specific primers. Approximately, one microgram of RNA is reverse transcribed with MolMRT (reverse transcriptase) and DDB (CDHR20), DQB (CDQB07), DQA (CDQA05), DPB (DPB11) and DPA (DPA14, DPA13) (optional) -specific non-sense primers in a 20 μ l final volume reaction (30-60 minute incubation). The

reaction for each Class II gene is performed in a different tube, but they can be performed in the same tube if performed. For routine purposes, four simultaneous reactions are performed for DQB, one for DQB, one for DQA, one for DPA1, and two for DPA1 gene products.

Once these reactions are completed, the enzymatic amplification of the respective cDNA molecules is then performed by directly adding to the 20 μ l reverse transcription reaction, the reagents needed for the amplification step. Alternatively, if DNA is used, the primer combinations used for the PCR are those shown in Table II herein (the anti-sense primers as well as the sense primers will be different). This includes the PCR reagents and appropriate conserved and non-conserved oligonucleotide primers. This example uses four reactions for DQB (tubes 1, 2, 3 and 4), one for DQB (tube 5), one for DQA (tube 6), one for DQB (tube 7), and two for DPA (tubes 8 and 9, respectively). Reactions 2, 3 and 4 incorporate primers DQB23, DQB24 and DQB25, respectively. For rapid typing (in less than 24 hours), the latter are the preferred combinations. Alternative combinations of the primers that can be used are shown in Table II.

Once completed, the reactions are spun-dialyzed for about 15 minutes using Centricon (Amicon, Ultraclean (mili-pore)) or similar columns to remove unincorporated primers and dNTPs. The retentate or one half of the recovered reactants for each reaction is then directly sequenced using Taq polymerase and the primers described in Table II for each combination of primers used in the cDNA/PCR reactions using 30 p-32 end-labeled (10 minutes) locus-specific sequencing primers (30 minutes).

The sequencing reactions products are loaded on an acrylamide gel, electrophoresed in 2-3 hours and exposed to X-ray films for 4-12 hours. The gels are read and results from

gels are compared to nucleotide sequences corresponding to all possible alleles.

Comparisons can be made visually using the naked eye or using a personal computer and a software package including the nucleotide sequences of all alleles of all haplotypes and routines which indicate how the comparison is to be performed as well as subroutines which will allow identification of new allelic sequences.

WO 92/15711

PCT/US92/01675

-52-

SEQUENCE LISTING

- (1) GENERAL INFORMATION:
- (i) APPLICANT: Santamaría, Pedro
Boyes-Jacino, Michael T.
Barbosa, Jose J.
Rich, Stephen S.
Zare, Anthony J.
- (ii) TITLE OF INVENTION: DNA Sequence-Based HLA Typing Method
- (iii) NUMBER OF SEQUENCES: 49
- (iv) CORRESPONDENCE ADDRESS:
- (A) ADDRESSEE: Merchant & Gould
(B) STREET: 3100 Northwest Center
(C) CITY: Minneapolis
(D) STATE: Minnesota
(E) COUNTRY: USA
(F) ZIP: 55402
- (v) COMPUTER READABLE FORM:
- (A) MEDIUM TYPE: Diskette, 3.5 inch,
720 Kb
(B) COMPUTER: Northgate 386
(C) OPERATING SYSTEM: DOS 4.0
(D) SOFTWARE: WordPerfect 5.0
- (vi) CURRENT APPLICATION DATA:
- (A) APPLICATION NUMBER: 07/665,960
(B) FILING DATE: 06-MAR-1991
(C) CLASSIFICATION:
- (vii) ATTORNEY INFORMATION:
- (A) NAME: Kowalsky, Allen W.
(B) REGISTRATION NUMBER: 31,435
(C) REFERENCE/DOCKET NUMBER: 600.24JUS
600.24JWO
- (x) TELECOMMUNICATION INFORMATION:
- (A) TELEPHONE: (612) 332-3300
(B) TELEFAX: (612) 332-9081

(2) INFORMATION FOR SEQUENCE ID NO: 1:

- (1) SEQUENCE CHARACTERISTICS:
 (A) LENGTH: 21 base pairs
 (B) TYPE: Nucleic Acid
 (C) STRANDEDNESS: Single
 (D) TOPOLOGY: Linear
- (11) MOLECULE TYPE: Genomic DNA
- (1v) ANTI-SENSE: Yes
- (v) FRAGMENT TYPE: Internal Fragment
- (v1) ORIGINAL SOURCE: Synthetically Derived
- (1x) FEATURE:
- (A) NAME/KEY: Oligonucleotide Primer DQB7
- (B) LOCATION: Anneals to codons 105 to 111 of the DQB1 transcript of HLA class II
- (x1) SEQUENCE DESCRIPTION: SEQ ID NO: 1:
 G GTC GTT GAG GGC CTC TGT CC 21
- (2) INFORMATION FOR SEQUENCE ID NO: 2:
- (1) SEQUENCE CHARACTERISTICS:
 (A) LENGTH: 21 base pairs
 (B) TYPE: Nucleic Acid
 (C) STRANDEDNESS: Single
 (D) TOPOLOGY: Linear
- (11) MOLECULE TYPE: Genomic DNA
- (1v) ANTI-SENSE: Yes
- (v) FRAGMENT TYPE: Internal Fragment
- (v1) ORIGINAL SOURCE: Synthetically Derived
- (1x) FEATURE:
- (A) NAME/KEY: Oligonucleotide Primer DRB20

- (B) LOCATION: Anneals to codons 105 to 111 of the DRB1, DRB3, DRB4 and DRB5 transcripts of HLA class II
- (x1) SEQUENCE DESCRIPTION: SEQ ID NO: 2:
 GAG CTG CAG GGG CTG GGT CTT 21
- (2) INFORMATION FOR SEQUENCE ID NO: 3:
- (1) SEQUENCE CHARACTERISTICS:
 (A) LENGTH: 22 base pairs
 (B) TYPE: Nucleic Acid
 (C) STRANDEDNESS: Single
 (D) TOPOLOGY: Linear
- (11) MOLECULE TYPE: Genomic DNA
- (1v) ANTI-SENSE: Yes
- (v) FRAGMENT TYPE: Internal Fragment
- (v1) ORIGINAL SOURCE: Synthetically Derived
- (1x) FEATURE:
- (A) NAME/KEY: Oligonucleotide Primer DQA9
- (B) LOCATION: Anneals to codons 148 to 155 of the DQA1 transcript of HLA class II
- (x1) SEQUENCE DESCRIPTION: SEQ ID NO: 3:
 GGT GAG GTT ACT GAT CTT GAA G 22
- (2) INFORMATION FOR SEQUENCE ID NO: 4:
- (1) SEQUENCE CHARACTERISTICS:
 (A) LENGTH: 21 base pairs
 (B) TYPE: Nucleic Acid
 (C) STRANDEDNESS: Single
 (D) TOPOLOGY: Linear

-55-

(11) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: no

(v) FRAGMENT TYPE: Internal Fragment

5 ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

10 (A) NAME/KEY: Oligonucleotide Primer:
DQB13

15 (B) LOCATION: Anneals to codons 1 to 7 of
the DQB1 transcript of HLA class II

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 4:

20 AGA GAC TCT CCC GAG GAT TTC 21
Arg Asp Ser Pro Glu Asp Phe
1 5

25 (2) INFORMATION FOR SEQUENCE ID NO: 5:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 21 base pairs

(B) TYPE: Nucleic Acid

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

30 (A) NAME/KEY: Oligonucleotide Primer:
DRB22

35 (A) LENGTH: 21 base pairs

(B) TYPE: Nucleic Acid

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: no

(v) FRAGMENT TYPE: Internal Fragment

40 (vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

45 (A) NAME/KEY: Oligonucleotide Primer:
DRB22

50 (B) LOCATION: Anneals to codons -4 to +3
of the DRB1, DRB3, DRB4 and DRB5
transcripts of HLA class II

-56-

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 5:

CTG GCT TTG GGT GGG GAC ACC 21
Leu Ala Leu Ala Gly Asp Thr
1

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 21 base pairs

(B) TYPE: Nucleic Acid

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

15 (ii) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: no

(v) FRAGMENT TYPE: Internal Fragment

20 (vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

25 (A) NAME/KEY: Oligonucleotide Primer:
DRB11

30 (B) LOCATION: Anneals to codons -33 to -
26 of the DRB1, DRB3, DRB4 and DRB5
transcripts of HLA class II

35 (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 6:

40 TGT TTC ACC ACC ATG GAG TGT C 21
Phe Ser Ser Ser Val Cys Leu
-30

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 20 base pairs

(B) TYPE: Nucleic Acid

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

45 (ii) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: no

-59-

(ix)

FEATURE:

(A) NAME/KEY: Oligonucleotide Primer
DQB30

5

(B) LOCATION: Anneals to codons 97 to 104
of the DQB1 transcript of HLA class
II

10 (xi)

SEQUENCE DESCRIPTION: SEQ ID NO: 10:

A TCG GGA GAT GGT CAC TGT GG 21

15

(2) INFORMATION FOR SEQUENCE ID NO: 11:

(i)

SEQUENCE CHARACTERISTICS:

(A) LENGTH: 21 base pairs

(B) TYPE: Nucleic Acid

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

25

(ii) MOLECULE TYPE: Genomic DNA

(iv)

ANTI-SENSE: Yes

30 (v)

FRAGMENT TYPE: Internal Fragment

(vi) ORIGINAL SOURCE: Synthetically Derived

35 (ix)

FEATURE:

(A) NAME/KEY: Oligonucleotide Primer
DRB30

40

(B) LOCATION: Anneals to codons 97 to 103
of the DRB1, DRB3, DRB4 and DRB5
transcripts of HLA class II

45 (xi)

SEQUENCE DESCRIPTION: SEQ ID NO: 11:

AGG ATA CAC AGT CAC CTT AGG 21

50 (2)

INFORMATION FOR SEQUENCE ID NO: 12:

(i)

SEQUENCE CHARACTERISTICS:

(A) LENGTH: 18 base pairs

(B) TYPE: Nucleic Acid

55

WO 92/15711

PCT/US92/01675

-60-

(C) STRANDEDNESS: Single
(D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: Yes

(v) FRAGMENT TYPE: Internal Fragment

(vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

(A) NAME/KEY: Oligonucleotide Primer DQB5

(B) LOCATION: Anneals to codons 78 to 83
of the DQB1 transcript of HLA class
II

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 12:

GTA GTT GTG TCT GCA CAC 18

(2) INFORMATION FOR SEQUENCE ID NO: 13:

(i)

SEQUENCE CHARACTERISTICS:

(A) LENGTH: 20 base pairs

(B) TYPE: Nucleic Acid

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

(ii)

MOLECULE TYPE: Genomic DNA

(iv)

ANTI-SENSE: Yes

(v)

FRAGMENT TYPE: Internal Fragment

(vi)

ORIGINAL SOURCE: Synthetically Derived

(ix)

FEATURE:

(A) NAME/KEY: Oligonucleotide Primer
DRB12(B) LOCATION: Anneals to codons 87 to 94
of the DRB1, DRB3, DRB4 and DRB5
transcripts of HLA class II

WO 92/15711

PCT/US92/01675

-61-

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 13:

G CCG CTG CAC TGT GAA GGT C 20

5

(2) INFORMATION FOR SEQUENCE ID NO: 14:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 23 base pairs

(B) TYPE: Nucleic Acid

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: yes

(v) FRAGMENT TYPE: Internal Fragment

(vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

(A) NAME/KEY: Oligonucleotide Primer

DQA25

(B) LOCATION: Anneals to codons 82 to 89
of the DQA1 transcript of HLA class II

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 14:

CAC GGT TCC GGT ACC ACC GGT AG 23

40

(2) INFORMATION FOR SEQUENCE ID NO: 15:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 18 base pairs

(B) TYPE: Nucleic Acid

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: no

(v) FRAGMENT TYPE: Internal Fragment

55

WO 92/15711

PCT/US92/01675

-62-

(vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

(A) NAME/KEY: Oligonucleotide Primer
DQA30(B) LOCATION: Anneals to codons 19 to 24
of the DQA1 transcript of HLA class II

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 15:

TAC GGT CCC TGT GTC CAG 18
ATG GAT TTT GGC GAT GAT 20

(2) INFORMATION FOR SEQUENCE ID NO: 16:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 20 base pairs

(B) TYPE: Nucleic Acid

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: no

(v) FRAGMENT TYPE: Internal Fragment

(vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

(A) NAME/KEY: Oligonucleotide Primer
DRE1400(B) LOCATION: Anneals to codons 38 to 45
of the DRE1 transcript of HLA class II

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 16:

G CGC TTC GAC ACC GAC GTC G 20
Val Arg Phe Asp Ser Asp Val Gly 40 45

(2) INFORMATION FOR SEQUENCE ID NO: 17:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 21 base pairs
(B) TYPE: Nucleic Acid
(C) STRANDEDNESS: Single
(D) TOPOLOGY: Linear

5 (ii) MOLECULE TYPE: Genomic DNA

10 (iv) ANTI-SENSE: no

(v) FRAGMENT TYPE: Internal Fragment

15 (vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

20 (A) NAME/KEY: Oligonucleotide Primer
DBR1401

(B) LOCATION: Anneals to codons 98 to 104
of the DBR1-0701-2 transcript of H1A
class II

25 (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 17:

GAG GTG ACT GTG TAT CCT GAC 21
Glu Val Thr Val Tyr Pro Asp
100

35 (2) INFORMATION FOR SEQUENCE ID NO: 18:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 21 base pairs
(B) TYPE: Nucleic Acid
(C) STRANDEDNESS: Single
(D) TOPOLOGY: Linear

40 (ii) MOLECULE TYPE: Genomic DNA

45 (iv) ANTI-SENSE: yes

(v) FRAGMENT TYPE: Internal Fragment

50 (vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

(A) NAME/KEY: Oligonucleotide Primer
DQB1402

(B) LOCATION: Anneals to codons 142 to
148 of the DBR1, DBR3, DBR4 and DBR5
transcripts of H1A class II

5 (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 18:

GAT CAG GGC TGT GCA CAC CAC 21

10 (2) INFORMATION FOR SEQUENCE ID NO: 19:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 21 base pairs
(B) TYPE: Nucleic Acid
(C) STRANDEDNESS: Single
(D) TOPOLOGY: Linear

20 (ii) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: yes

(v) FRAGMENT TYPE: Internal Fragment

(vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

(A) NAME/KEY: Oligonucleotide Primer
DQB1403

35 (B) LOCATION: Anneals to codons 127 to
133 of the DBR1, DBR3, DBR4 and DBR5
transcripts of H1A class II

40 (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 19:

CCG GAA CCA CCT GAC TTC AAT 21

45 (2) INFORMATION FOR SEQUENCE ID NO: 20:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 21 base pairs
(B) TYPE: Nucleic Acid
(C) STRANDEDNESS: Single
(D) TOPOLOGY: Linear

55 (ii) MOLECULE TYPE: Genomic DNA

- (iv) ANTI-SENSE: NO
- (v) FRAGMENT TYPE: Internal Fragment
- 5 (vi) ORIGINAL SOURCE: Synthetically Derived
- (ix) FEATURE:
- 10 (A) NAME/KEY: Oligonucleotide Primer
DOB1406
- (B) LOCATION: Anneals to bp18-38 to
int. 33 of the DBR1, DBR3, DBR4 and
DBR5 transcripts of HLA class II
- 15 (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 20:
GCCACGATG GGCCTGCGC C 21
- 20
- 25 (2) INFORMATION FOR SEQUENCE ID NO: 21:
- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 21 base pairs
- (B) TYPE: Nucleic Acid
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: Genomic DNA
- 30 (iv) ANTI-SENSE: Yes
- (v) FRAGMENT TYPE: Internal Fragment
- (vi) ORIGINAL SOURCE: Synthetically Derived
- (ix) FEATURE:
- 40 (A) NAME/KEY: Oligonucleotide Primer
DBR825
- (B) LOCATION: Anneals to codons 79 to 85
of the DBR1, DBR3, DBR4 and DBR5
transcripts of HLA class II
- 45
- 50 (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 21:
AAC CCC GTA GTT GTG TCT GCA 21
- 55

- (2) INFORMATION FOR SEQUENCE ID NO: 22:
- 5 (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 21 base pairs
- (B) TYPE: Nucleic Acid
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: Genomic DNA
- (iv) ANTI-SENSE: NO
- (v) FRAGMENT TYPE: Internal Fragment
- 10 (vi) ORIGINAL SOURCE: Synthetically Derived
- (ix) FEATURE:
- (A) NAME/KEY: Oligonucleotide Primer
DBR824
- (B) LOCATION: Anneals to codons 1 to 7 of
the DBR1, DBR3, DBR4 and DBR5
transcripts of HLA class II
- 25 (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 22:
GGG GAC ACC CGA CCA CGT TTC 21
1
Cly Ala Thr Arg Pro Arg Phe
- 30
- 35
- (2) INFORMATION FOR SEQUENCE ID NO: 23:
- 40 (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 21 base pairs
- (B) TYPE: Nucleic Acid
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: Genomic DNA
- (iv) ANTI-SENSE: NO
- (v) FRAGMENT TYPE: Internal Fragment
- 45
- 50 (vi) ORIGINAL SOURCE: Synthetically Derived
- (ix) FEATURE:
- 55

-67-

- (A) NAME/KEY: Oligonucleotide Primer
DPB10
- (B) LOCATION: Anneals to codons -19 to -13 of the DPB1 transcript of HLA class II
- (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 23:
CGG ACA GTG GCT CTG ACG GCG 21
Arg Thr Val Ala Leu Tyr Ala
-15
- (2) INFORMATION FOR SEQUENCE ID NO: 24:
- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 21 base pairs
- (B) TYPE: Nucleic Acid
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: Genomic DNA
- (iv) ANTI-SENSE: yes
- (v) FRAGMENT TYPE: Internal Fragment
- (vi) ORIGINAL SOURCE: Synthetically Derived
- (ix) FEATURE:
- (A) NAME/KEY: Oligonucleotide Primer
DPB11
- (B) LOCATION: Anneals to codons 105 to 111 of the DPB1 transcript of HLA class II
- (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 24:
CTT GTG GTG CTG CAA GGG CCC 21
- (2) INFORMATION FOR SEQUENCE ID NO: 25:
- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 21 base pairs
- (B) TYPE: Nucleic Acid
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: Genomic DNA
- (iv) ANTI-SENSE: yes
- (v) FRAGMENT TYPE: Internal Fragment
- (vi) ORIGINAL SOURCE: Synthetically Derived
- (ix) FEATURE:
- (A) NAME/KEY: Oligonucleotide Primer
DPB13
- (B) LOCATION: Anneals to codons -5 to +2 of the DPB1 transcript of HLA class II

-68-

- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: Genomic DNA
- (iv) ANTI-SENSE: yes
- (v) FRAGMENT TYPE: Internal Fragment
- (vi) ORIGINAL SOURCE: Synthetically Derived
- (ix) FEATURE:
- (A) NAME/KEY: Oligonucleotide Primer
DPB12
- (B) LOCATION: Anneals to codons 97 to 103 of the DPB1 transcript of HLA class II
- (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 25:
CTT GGA GCG GCA AAC ATT CAC 21
- (2) INFORMATION FOR SEQUENCE ID NO: 26:
- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 21 base pairs
- (B) TYPE: Nucleic Acid
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: Genomic DNA
- (iv) ANTI-SENSE: no
- (v) FRAGMENT TYPE: Internal Fragment
- (vi) ORIGINAL SOURCE: Synthetically Derived
- (ix) FEATURE:
- (A) NAME/KEY: Oligonucleotide Primer
DPB13
- (B) LOCATION: Anneals to codons -5 to +2 of the DPB1 transcript of HLA class II

-69-

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 26:

TA CTC AGC GTC CTC CTC ACA T

Lau Lau Met Val Lau Lau Thr Ser

-12 -5

5

(2) INFORMATION FOR SEQUENCE ID NO: 27:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 21 base pairs

(B) TYPE: Nucleic Acid

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: no

(v) FRAGMENT TYPE: Internal Fragment

(vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

(A) NAME/KEY: Oligonucleotide Primer
DPB14

30

(B) LOCATION: Anneals to bp-42/-46 to
intron 2 of the DPB1 transcript of
HLA class II

SEQUENCE DESCRIPTION: SEQ ID NO: 27:

AGAGGAGAA AGAGTATG A

21

35

(2) INFORMATION FOR SEQUENCE ID NO: 28:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 21 base pairs

(B) TYPE: Nucleic Acid

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: no

(v) FRAGMENT TYPE: Internal Fragment

(vi) ORIGINAL SOURCE: Synthetically Derived

55

-70-

(ix)

FEATURE:

(A) NAME/KEY: Oligonucleotide Primer
DPB15

5

(B) LOCATION: Anneals to bp39-59 to
intron 3 of the DPB1 transcript of
HLA class II

SEQUENCE DESCRIPTION: SEQ ID NO: 28:

GCCCTGAGCA CGGGCCCGG C

21

15

(2) INFORMATION FOR SEQUENCE ID NO: 29:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 21 base pairs

(B) TYPE: Nucleic Acid

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: no

(v) FRAGMENT TYPE: Internal Fragment

(vi) ORIGINAL SOURCE: Synthetically Derived

(ix)

FEATURE:

(A) NAME/KEY: Oligonucleotide Primer
DPB16

40

(B) LOCATION: Anneals to bp1-21 to intron
3 of the DPB1 transcript of HLA class
II

SEQUENCE DESCRIPTION: SEQ ID NO: 29:

CGGCCCAAG CCTGACCTA C

21

45

50

(2) INFORMATION FOR SEQUENCE ID NO: 30:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 21 base pairs

55

-71-

(B) TYPE: Nucleic Acid
(C) STRANDEDNESS: Single
(D) TOPOLOGY: Linear

(11) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: no

(v) FRAGMENT TYPE: Internal Fragment

(vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

(A) NAME/KEY: Oligonucleotide Primer
DP917

(B) LOCATION: Anneals to bp-5/-26 to
110 of the DP91 transcript of
HLA class II

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 30:

CGCTCAAGTC CCCCCCTCC C 21

(2) INFORMATION FOR SEQUENCE ID NO: 31:

(1) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 21 base pairs

(B) TYPE: Nucleic Acid

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

(11) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: yes

(v) FRAGMENT TYPE: Internal Fragment

(vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

(A) NAME/KEY: Oligonucleotide Primer
DP914

(B) LOCATION: Anneals to codons 104 to
110 of the DP91 transcript of HLA
class II

55

-72-

SEQUENCE DESCRIPTION: SEQ ID NO: 31:

GTG AAT GTG GCA GAT GAG GGT 21

5

(2) INFORMATION FOR SEQUENCE ID NO: 32:

(1) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 21 base pairs

(B) TYPE: Nucleic Acid

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

(11) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: yes

(v) FRAGMENT TYPE: Internal Fragment

(vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

(A) NAME/KEY: Oligonucleotide Primer
DP915

(B) LOCATION: Anneals to codons -17 to -
23 of the DP91 transcript of HLA
class II

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 32:

CAT ATC AGA GGT GTG ATC TTG 21

(2) INFORMATION FOR SEQUENCE ID NO: 33:

(1) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 21 base pairs

(B) TYPE: Nucleic Acid

(C) STRANDEDNESS: Single

(D) TOPOLOGY: Linear

(11) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: yes

(v) FRAGMENT TYPE: Internal Fragment

-73-

- (vi) ORIGINAL SOURCE: Synthetically Derived
- (ix) FEATURE:
- (A) NAME/KEY: Oligonucleotide Primer
DPA16
- (B) LOCATION: Anneals to codons 88 to 94
of the DPA1 transcript of H1A class II
- (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 33:
- CTT GGG AAA CAC GGT CAC CTC 21
- (2) INFORMATION FOR SEQUENCE ID NO: 34:
- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 21 base pairs
- (B) TYPE: Nucleic Acid
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: Genomic DNA
- (iv) ANTI-SENSE: no
- (v) FRAGMENT TYPE: Internal Fragment
- (vi) ORIGINAL SOURCE: Synthetically Derived
- (ix) FEATURE:
- (A) NAME/KEY: Oligonucleotide Primer
DPA17
- (B) LOCATION: Anneals to codons -3 to -9
of the DPA1 transcript of H1A class II
- (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 34:
- CTG CTG ACT CTC GGA GGA GCT 21
Leu Leu Ser Leu Arg Gly Ala
-5

55

-74-

- (2) INFORMATION FOR SEQUENCE ID NO: 35:
- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 20 base pairs
- (B) TYPE: Nucleic Acid
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: Genomic DNA
- (iv) ANTI-SENSE: no
- (v) FRAGMENT TYPE: Internal Fragment
- (vi) ORIGINAL SOURCE: Synthetically Derived
- (ix) FEATURE:
- (A) NAME/KEY: Oligonucleotide Primer
DPA18
- (B) LOCATION: Anneals to bp-69/-50 of
intron 2 of the DPA1 transcript of
H1A class II
- (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 35:
- CTCTACCTT GACCACTTC 20
- (2) INFORMATION FOR SEQUENCE ID NO: 36:
- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 21 base pairs
- (B) TYPE: Nucleic Acid
- (C) STRANDEDNESS: Single
- (D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: Genomic DNA
- (iv) ANTI-SENSE: no
- (v) FRAGMENT TYPE: Internal Fragment
- (vi) ORIGINAL SOURCE: Synthetically Derived
- (ix) FEATURE:
- (A) NAME/KEY: Oligonucleotide Primer
DPA19

55

-75-

(B) LOCATION: Anneals to bp55-71 to Intron 3 of the DPAl transcript of HLA class II

5 (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 36: 21
ACTCTGAGG TGGCAGAG G

10 (2) INFORMATION FOR SEQUENCE ID NO: 37:

(1) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 21 base pairs
(B) TYPE: Nucleic Acid
(C) STRANDEDNESS: Single
(D) TOPOLOGY: Linear

20 (ii) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: yes

(v) FRAGMENT TYPE: Internal Fragment

(vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

(A) NAME/KEY: Oligonucleotide Primer
DPAl8

35 (B) LOCATION: Anneals to codons 76 to 82 of the DPAl transcript of HLA class II

40 (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 37: 21
GGC CTG AGT GTG GTT GGA ACG

45 (2) INFORMATION FOR SEQUENCE ID NO: 38:

(1) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 21 base pairs
(B) TYPE: Nucleic Acid
(C) STRANDEDNESS: Single
(D) TOPOLOGY: Linear

50 (ii) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: no

-76-

(v) FRAGMENT TYPE: Internal Fragment

(vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

(A) NAME/KEY: Oligonucleotide Primer
DPAl8

(B) LOCATION: Anneals to codons 59 to 65 of the DPAl transcript of HLA class II

15 (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 38:

CTG GCT AAC ATT GCT AFA TGG 21
Leu Ala Asn Ile Ala Ile Leu 65

25 (2) INFORMATION FOR SEQUENCE ID NO: 39:

(1) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 21 base pairs
(B) TYPE: Nucleic Acid
(C) STRANDEDNESS: Single
(D) TOPOLOGY: Linear

(ii) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: yes

(v) FRAGMENT TYPE: Internal Fragment

(vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

(A) NAME/KEY: Oligonucleotide Primer
DPAl9

(B) LOCATION: Anneals to codons 222 to 228 of the DPAl transcript of HLA class II

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 39:

GGT CCC CTG GGC CCC GGG GTC 21

-77-

- 5 (2) INFORMATION FOR SEQUENCE ID NO. 40:
- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 21 base pairs
(B) STRANDS: 2
(C) STRANDEDNESS: Single
(D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: Genomic DNA
- (iv) ANTI-SENSE: yes
- (v) FRAGMENT TYPE: Internal Fragment
- (vi) ORIGINAL SOURCE: Synthetically Derived
- (ix) FEATURE:
- (A) NAME/KEY: Oligonucleotide Primer
DPA20
- (B) LOCATION: Anneals to codons 214 to 220 of the DPAL transcript of HLA class II
- (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 40:
- GCC AGA ACG CAG AGA CTT TAT 21
- 35 (2) INFORMATION FOR SEQUENCE ID NO. 41:
- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 21 base pairs
(B) STRANDS: 2
(C) STRANDEDNESS: Single
(D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: Genomic DNA
- (iv) ANTI-SENSE: no
- (v) FRAGMENT TYPE: Internal Fragment
- (vi) ORIGINAL SOURCE: Synthetically Derived

-78-

- (ix) FEATURE:
- (A) NAME/KEY: Oligonucleotide Primer
DPA21
- (B) LOCATION: Anneals to codons 66 to 74 of the DPAL transcript of HLA class II
- (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 41:
- AAC TTG AAT ACC TTG ATC CAG 21
Asn Leu Asn Thr Leu Ile Gln 70
- 15 (2) INFORMATION FOR SEQUENCE ID NO. 42:
- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 20 base pairs
(B) STRANDS: 2
(C) STRANDEDNESS: Single
(D) TOPOLOGY: Linear
- (ii) MOLECULE TYPE: Genomic DNA
- (iv) ANTI-SENSE: no
- (v) FRAGMENT TYPE: Internal Fragment
- (vi) ORIGINAL SOURCE: Synthetically Derived
- (ix) FEATURE:
- (A) NAME/KEY: Oligonucleotide Primer
DPA23
- (B) LOCATION: Anneals to codons 7 to 13 of the DRB1 transcript of HLA class II
- (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 42:
- ATC TTG CAG CAG CAT AAG TA 20
Phe Leu Gln Cys Asp Lys Tyr 10
- 50

-79-

(2) INFORMATION FOR SEQUENCE ID NO: 43:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 21 base pairs
 (B) TYPE: Nucleic Acid
 (C) STRANDEDNESS: Single
 (D) TOPOLOGY: Linear

5

(ii) MOLECULE TYPE: Genomic DNA

10

(iv) ANTI-SENSE: no

15

(v) FRAGMENT TYPE: Internal Fragment

(vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

- (A) NAME/KEY: Oligonucleotide Primer
 DR824

20

- (B) LOCATION: Annals to codons 5 to 11
 of the DRB1 transcript of HLA class
 II

25

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 43:

CCA CGT TTC TGC GAC TAC TCT 21
 Pro Arg Phe Leu Gly Tyr Ser
 5

30

(2) INFORMATION FOR SEQUENCE ID NO: 44:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 21 base pairs
 (B) TYPE: Nucleic Acid
 (C) STRANDEDNESS: Single
 (D) TOPOLOGY: Linear

40

(ii) MOLECULE TYPE: Genomic DNA

45

(iv) ANTI-SENSE: no

50

(v) FRAGMENT TYPE: Internal Fragment

(vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

- (A) NAME/KEY: Oligonucleotide Primer
 DR825

55

-80-

- (B) LOCATION: Annals to codons 6 to 13
 of the DRB1 transcript of HLA class
 II

5

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 44:

T TTC TTC GAG CAG GGT AAA CA 21
 Arg Phe Leu Glu Gln Val Lys His
 10

10

(2) INFORMATION FOR SEQUENCE ID NO: 45:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 21 base pairs
 (B) TYPE: Nucleic Acid
 (C) STRANDEDNESS: Single
 (D) TOPOLOGY: Linear

15

(ii) MOLECULE TYPE: Genomic DNA

25

(iv) ANTI-SENSE: no

30

(v) FRAGMENT TYPE: Internal Fragment

(vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

- (A) NAME/KEY: Oligonucleotide Primer
 DR816

35

- (B) LOCATION: Annals to codons 29 to 35
 of the DRB1 transcript of HLA class II

40

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 45:

AGA TGC ATC TAT AAC CAA GAG 21
 Arg Cys Ile Tyr Asn Glu Gln
 30

45

(2) INFORMATION FOR SEQUENCE ID NO: 46:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 21 base pairs
 (B) TYPE: Nucleic Acid

50

55

-81-

(C) STANDARDS: Single
(D) TOPOLOGY: Linear

(11) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: no

(v) FRAGMENT TYPE: Internal Fragment

(vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

(A) NAME/KEY: Oligonucleotide Primer
DBA17

(B) LOCATION: Anneals to codons 29 to 35
of the DBA1 DBA3 DBA5 DBA5
transcripts of HLA class II

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 46:

AGA AAC TTC CAT AAC CAG GAG 21
ATG TTT Phe His Asn Gln Glu 35
30

30 (2) INFORMATION FOR SEQUENCE ID NO: 47:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 21 base pairs

(B) TYPE: Nucleic Acid

(C) STANDARDS: Single

(D) TOPOLOGY: Linear

(11) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: no

(v) FRAGMENT TYPE: Internal Fragment

(vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

(A) NAME/KEY: Oligonucleotide Primer DBA6

(B) LOCATION: Anneals to codons -8 to -2
of the DBA1 transcript of HLA class
II

55

-82-

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 47:

CTG AGC ACC CCA GTG GCT GAG 21
Leu Ser Thr Pro Val Ala Glu
-5

5

10 (2) INFORMATION FOR SEQUENCE ID NO: 48:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 21 base pairs

(B) TYPE: Nucleic Acid

(C) STANDARDS: Single

(D) TOPOLOGY: Linear

(11) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: no

(v) FRAGMENT TYPE: Internal Fragment

(vi) ORIGINAL SOURCE: Synthetically Derived

(ix) FEATURE:

(A) NAME/KEY: Oligonucleotide Primer
DBA14

(B) LOCATION: Anneals to codons -8 to -2
of the DBA1 transcript of HLA class
II

(xi) SEQUENCE DESCRIPTION: SEQ ID NO: 48:

CTG ACC TCC TCA CTG GCT GAG 21
Leu Ser Ser Ser Leu Ala Glu
-5

45 (2) INFORMATION FOR SEQUENCE ID NO: 49:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 21 base pairs

(B) TYPE: Nucleic Acid

(C) STANDARDS: Single

(D) TOPOLOGY: Linear

(11) MOLECULE TYPE: Genomic DNA

(iv) ANTI-SENSE: no

-83-

- (v) FRAGMENT TYPE: Internal Fragment
- (vi) ORIGINAL SOURCE: Synthetically Derived
- (ix) FEATURE:
- (A) NAME/KEY: Oligonucleotide Primer
DDB15
- (B) LOCATION: Anneals to codons -8 to -2
in the Dgat transcript of Rha class
II
- (xi) SEQUENCE DESCRIPTION: SEQ ID NO: 49:
- CTG AGC AGC TCG CTG GCT GAG 21
Leu Ser Thr Ser Val Ala Glu
-5

Applicants state that the paper copy of the above "Sequences Listing" Section of the present application, and the computer readable form of the same submitted therewith, are the same.

84

WHAT IS CLAIMED:

1. A method for determining a major histocompatibility complex genotype of a subject in a sample containing subject nucleic acid comprising:
 - (a) isolating nucleic acid from said sample;
 - (b) amplifying said nucleic acid by polymerase chain reaction to generate sufficient polymerase chain reaction product for each allele of said gene locus to be sequenced, all or said alleles for each gene locus and chromosome to be sequenced being amplified with at least one conserved oligonucleotide primer pair and at least one of said alleles for each gene locus and chromosome to be sequenced being amplified with at least one conserved oligonucleotide primer and at least one non-conserved oligonucleotide primer;
 - (c) sequencing directly each polymerase chain reaction product for each allele at each gene locus of each chromosome with Tag polymerase and a conserved primer specific for each locus that is sequenced; and
 - (d) analyzing each sequenced polymerase chain reaction product to determine the genotype of said subject.
2. The method of claim 1 wherein said isolated nucleic acid is genomic DNA.
3. The method of claim 1 wherein said isolated nucleic acid is RNA and further comprises the following step prior to amplifying said nucleic acid:
 - (a) synthesizing cDNA molecules for each allele of each gene locus to be sequenced, wherein said

synthesis employs a locus-specific oligonucleotide primer that anneals to a conserved region of each allele of each said gene locus.

4. The method of claim 1 wherein said major histocompatibility genotype to be determined is a HLA Class II genotype.

5. The method of claim 4 wherein said Class II gene locus to be sequenced is DQB1.

6. The method of claim 4 wherein said Class II gene locus to be sequenced is DQA1.

7. The method of claim 4 wherein said Class II gene loci to be sequenced are DQB 1/3/4/5.

8. The method of claim 4 wherein said Class II gene loci to be sequenced is DPA 1.

9. The method of claim 4 wherein said Class II gene loci to be sequenced is DPA 1.

10. The method of claim 1 wherein analyzing said sequenced polymerase chain reaction product involves comparing the nucleotide sequence of each allele of each gene locus sequenced to known sequences for each such gene locus followed by comparing the sequence of each allele of each gene locus amplified with a conserved/non-conserved oligonucleotide primer pair to the nucleotide sequence of each allele of such gene locus amplified with a conserved oligonucleotide primer pair.

11. The method of claim 1 wherein analyzing each polymerase chain reaction product to determine genotype is conducted with a computer having a program including nucleotide sequences of all alleles of all haplotypes for HLA Class II loci.

12. The method of claim 1 wherein said amplifying cDNA molecules with said conserved oligonucleotide primer includes annealing said conserved oligonucleotide primer to said cDNA at about 37°C.

13. The method of claim 1 wherein said amplifying cDNA molecules with said non-conserved primer includes annealing said non-conserved primer to said cDNA at about 55°C.

14. A method for determining the Class II histocompatibility genotype of a subject in a sample containing subject nucleic acid comprising:

- (a) isolating nucleic acid from said sample;
- (b) amplifying said nucleic acid by polymerase chain reaction to generate sufficient polymerase chain reaction product for each allele of said Class II gene locus to be sequenced, all of said alleles for each Class II gene locus and chromosome to be sequenced being amplified with at least one conserved oligonucleotide primer pair and at least one of said alleles for each Class II gene locus and chromosome to be sequenced being amplified with at least one conserved oligonucleotide primer and at least one non-conserved oligonucleotide primer;
- (c) sequencing directly each polymerase chain reaction product for each allele at each Class

87

II gene locus of each chromosome with Tag polymerase and a conserved primer specific for each Class II locus that is sequenced; and

(d) determining the genotype of said subject by comparing the nucleotide sequence of each allele at each Class II locus sequenced to known sequences for each such Class II locus followed by comparing the sequence of each allele of each Class II locus amplified with a degenerated oligonucleotide primer to the nucleotide sequence of each allele of such Class II locus amplified with a conserved oligonucleotide primer.

15. The method of claim 14 wherein said isolated nucleic acid is RNA and further comprises the following step prior to amplifying said nucleic acid:

(a) synthesizing cDNA molecules for each allele of each Class II gene locus to be sequenced, wherein said synthesis employs a locus-specific oligonucleotide primer that anneals to a conserved region of each allele of each said Class II gene locus.

16. A method for determining the Class II HLA genotype of a subject in a sample containing subject nucleic acid comprising:

(a) isolating total cellular RNA from said sample; and

(b) synthesizing cDNA molecules for each allele of at least one Class II gene locus to be sequenced, wherein said synthesis employs a locus-specific oligonucleotide primer that anneals to a conserved region of each allele of each said Class II gene locus;

88

(c) amplifying said cDNA molecules by polymerase chain reaction to generate a polymerase chain reaction product for each allele of said Class II gene locus to be sequenced, all of said alleles for each Class II gene locus and chromosome to be sequenced being amplified with at least one conserved oligonucleotide primer pair and at least one of said alleles for each Class II gene locus and chromosome to be sequenced being amplified with at least one conserved oligonucleotide primer and at least one non-conserved oligonucleotide primer; sequencing directly each polymerase chain reaction product for each allele at each Class II gene locus of each chromosome with Tag polymerase and a conserved primer specific for each Class II locus that is sequenced to produce a nucleic acid sequence ladder for each allele; and

(d) analyzing each nucleic acid ladder to determine the genotype of said subject by comparing the nucleotide sequence of each allele of each Class II locus sequenced to known sequences for each such Class II locus followed by comparing the sequence of each allele of each Class II locus amplified with a conserved/non-conserved oligonucleotide primer pair to the nucleotide sequence of each allele of such Class II locus amplified with a conserved oligonucleotide primer pair.

17. The method of claim 16 wherein said Class II HLA genotype to be determined includes nucleotide sequences for the DRB1, DRB3, DRB4, DRB5, DQB1, DQA1 and DPB1 genes of said subject.

18. An oligonucleotide primer comprising a single strand of DNA which anneals to codons 105 to 111 of the DQB transcript.

19. An oligonucleotide primer having the sequence
GATGATGAGAGGCTCTGTGTC. (SEQ. ID NO:1)

20. An oligonucleotide primer comprising a single strand of DNA which anneals to codons 1 to 7 of DQB.

21. An oligonucleotide primer having the sequence
AGAGACTTCTCCGAGCTTTC. (SEQ. ID NO:4)

22. An oligonucleotide primer comprising a single strand of DNA which anneals to codons 148 to 155 of its DQA transcript.

23. An oligonucleotide primer having the sequence
GATGAGTCTACTGATCTTGAG. (SEQ. ID NO:13)

24. An oligonucleotide primer comprising a single strand of DNA which anneals to codons -10 to -4 of DQA cDNA.

25. An oligonucleotide primer having the sequence
CTGTCTCTCCGTGAGGCC. (SEQ. ID NO:17)

26. An oligonucleotide primer comprising a single strand of DNA which anneals to codons 105 to 111 of DBB1 transcript.

27. An oligonucleotide primer having the sequence
GTGCTGACAGGCTGTGGTCTT. (SEQ. ID NO:12)

28. An oligonucleotide primer comprising a single strand of DNA which anneals to codons -33 to -26 of DBB1 transcripts.

29. An oligonucleotide primer having the sequence
TCTTCTCAGAGAGCTGTGTC. (SEQ. ID NO:16)

30. An oligonucleotide primer comprising a single strand of DNA which anneals to codons 7 to 13 of DBB1 transcript.

31. An oligonucleotide primer having the sequence
TTCTTGACAGGATAGTA. (SEQ. ID NO:42)

32. An oligonucleotide primer comprising a single strand of DNA which anneals to codons 5 to 11 of DBB 1 transcript.

33. An oligonucleotide primer having the sequence
CCAGCTTCTTGAGTACTT. (SEQ. ID NO:43)

34. An oligonucleotide primer comprising a single strand of DNA which anneals to codons 6 to 13 of DBB 1 transcript.

35. An oligonucleotide primer having the sequence
TTCTTGAGAGGTTAACA. (SEQ. ID NO:44)

36. An oligonucleotide primer comprising a single strand of DNA which anneals to codons 105 to 111 of DBE transcript.

37. An oligonucleotide primer having the sequence
GTGAGGCTCTGAGGCC. (SEQ. ID NO:24)

39. An oligonucleotide primer comprising a single strand of DNA which anneals to codons 222 to 228 of DPA transcript.

39. An oligonucleotide primer having the sequence
GCTCCCTCCGCGCCCGGGGCTC. (SEQ. ID NO:39)

40. An oligonucleotide primer comprising a single strand of DNA which anneals to codons -19 to -13 of DPB transcript.

41. An oligonucleotide primer having the sequence
CTGACATGCTGCTCTGACGGCG. (SEQ. ID NO:23)

42. An oligonucleotide primer comprising a single strand of DNA which anneals to codons -23 to -17 of DPA transcript.

43. An oligonucleotide primer sequence of
CATATCTAGACGTGATCTTG. (SEQ. ID NO:32)

44. An oligonucleotide primer comprising a single strand of DNA which anneals to codons 59 to 65 of DPA transcript.

45. An oligonucleotide primer having the sequence
CTGCTTAACATGCTGATATTTG. (SEQ. ID NO:38)

46. An oligonucleotide primer comprising a single strand of DNA which anneals to codons 104 to 110 of DPA transcript.

47. An oligonucleotide primer having the sequence
GTCAATGTCGAGATGAGGT. (SEQ. ID NO:31)

48. An oligonucleotide primer having the sequence
GCCGCTGACATGGAAGCTC. (SEQ. ID NO:13)

49. An oligonucleotide primer having the sequence
CTGAGAGGAGGAAACATTCAC. (SEQ. ID NO:11)

50. An oligonucleotide primer having the sequence
GTATGTGTGTCGACAC. (SEQ. ID NO:12)

51. An oligonucleotide primer having the sequence
CACGGTTCGCGTAGCGCGTAG. (SEQ. ID NO:14)

52. An oligonucleotide primer having the sequence
CTTGAGGGGGAAACATTCAC. (SEQ. ID NO:25)

53. An oligonucleotide primer having the sequence
CTTGGGAAAACACGGTCACCTC. (SEQ. ID NO:33)

54. An oligonucleotide primer having the sequence
GCCGAGACGAGAGACTTAT. (SEQ. ID NO:40)

55. An oligonucleotide primer comprising a single strand of DNA which anneals to base pairs 18 to 38 of intron 3 of DEB loci.

56. An oligonucleotide primer having the sequence
of GTCAGACGTGGGCTTCGACC. (SEQ. ID NO:20)

57. An oligonucleotide primer comprising a single strand of DNA which anneals to codons -4 to +3 of the DEB transcript.

58. An oligonucleotide primer having the sequence
CTGGCTTGTGCTGGGACACC. (SEQ. ID NO:15)

59. An oligonucleotide primer comprising a single strand of DNA which anneals to codons 88 to 94 of the DQB transcript.

60. An oligonucleotide primer having the sequence
TCCCTCTCTCCAGGGTCCGGG. (SEQ. ID NO:8)

61. An oligonucleotide primer comprising a single strand of DNA which anneals to codons 11 to 17 of the DQB transcript.

62. An oligonucleotide primer having the sequence
of TTTAAGGCGCATGTCTACTTC. (SEQ. ID NO:9)

63. An oligonucleotide primer comprising a single strand of DNA which anneals to base pairs -42 to -62 of intron 2 of the DPB locus.

64. An oligonucleotide primer having the sequence
AGAGGAGGAAGAGCATTTAGA. (SEQ. ID NO:27)

65. An oligonucleotide primer comprising a single strand of DNA which anneals to intron 39 to 59 of intron 3 of the DPB gene.

66. An oligonucleotide primer having the sequence
GCCCTGGGACGGGCCGGG. (SEQ. ID NO:28)

67. An oligonucleotide primer comprising a single strand of DNA which anneals to base pairs -69 to -50 of intron 2 of the DPB locus.

68. An oligonucleotide primer having the sequence
CTCTACCTTCACCATTCG. (SEQ. ID NO:35)

69. An oligonucleotide primer comprising a single strand of DNA which anneals to base pairs 55 to 71 of intron 3 of the DPB locus.

70. An oligonucleotide primer having the sequence
AGCTTGAGGTGTCACAGAGG. (SEQ. ID NO:36)

71. An oligonucleotide primer comprising a single strand of DNA which anneals to codons 87 to 94 of the DRB transcript.

72. An oligonucleotide primer comprising a single strand of DNA which anneals to codons 38 to 45 of the DRB transcript.

73. An oligonucleotide primer having the sequence
GCCTCTGCACAGCAGCTGG. (SEQ. ID NO:16)

74. An oligonucleotide primer having the sequence
CGCCCTAAGGCTCTCTCTAC. (SEQ. ID NO:29)

75. An oligonucleotide primer having the sequence
GACCTGAGGTGGTGTGAAGG. (SEQ. ID NO:37)

76. An oligonucleotide primer having the sequence
TACTGATGGCTCTCTCTACAT. (SEQ. ID NO:26)

77. An oligonucleotide primer having the sequence
CGCTCAAGTCCGCCCCCTCC. (SEQ. ID NO:30)

78. An oligonucleotide primer having the sequence
CTCTGTAGTCTCCAGGAGCT. (SEQ. ID NO:34)

79. An oligonucleotide primer having the sequence
AACTGATATCTGTATCCG. (SEQ. ID NO:41)

80. An oligonucleotide primer having the sequence
ATGGGAGATGGTCACTGTGG. (SEQ. ID NO:10)

81. An oligonucleotide primer having the sequence
TACGCTCCTCTGGCAG. (SEQ. ID NO:15)

82. A method for rapid automated determination of
major histocompatibility complex class genotype of a subject
in a sample containing subject nucleic acid comprising:

- (a) isolating nucleic acid from said sample with an
RNA/DNA extractor;
- (b) amplifying said nucleic acid by polymerase
chain reaction using a thermocycler to generate
a polymerase chain reaction product for each
allele of each gene locus to be sequenced, all
of said alleles for each gene locus and
chromosome to be sequenced being amplified with
at least one conserved oligonucleotide primer
pair and at least one of said alleles for each
gene locus and chromosome to be sequenced being
amplified with at least one conserved
oligonucleotide primer and one non-conserved
oligonucleotide primer;
- (c) sequencing directly each polymerase chain
reaction product for each allele at each gene
locus of each chromosome in an automated
sequencing apparatus with Tag polymerase and a
conserved primer specific for each locus to be
sequenced; and

- (d) analyzing each sequenced polymerase chain
reaction product to determine the genotype of
said subject with a computer having a data base
with allelic sequence information to compare
the sequence of each allele of each gene locus

sequenced to known sequences for each such gene
locus followed by comparing the sequence of
each allele of each gene locus amplified with a
conserved/non-conserved oligonucleotide primer
pair to the nucleotide sequence of each allele
of such gene locus amplified with a conserved
oligonucleotide primer pair.

83. A method for determining the genotype at one or
more polymorphic gene locus of a subject in a sample
containing subject nucleic acid comprising:

- (a) isolating nucleic acid from said sample;
- (b) amplifying said nucleic acid by polymerase
chain reaction to generate sufficient
polymerase chain reaction product for each
allele of said gene locus to be sequenced, all
of said alleles for each gene locus and
chromosome to be sequenced being amplified with
at least one conserved oligonucleotide primer
pair and at least one of said alleles for each
gene locus and chromosome to be sequenced being
amplified with at least one conserved
oligonucleotide primer and at least one non-
conserved oligonucleotide primer;
- (c) sequencing directly each polymerase chain
reaction product for each allele at each gene
locus of each chromosome with a sequencing
enzyme and a conserved primer specific for each
locus that is sequenced; and
- (d) analyzing each sequenced polymerase chain
reaction product to determine the genotype of
said subject.

84. The method of claim 83 wherein said isolated
nucleic acid is genomic DNA.

97

85. The method of claim 83 wherein said isolated nucleic acid is RNA and further comprises the following step prior to amplifying said nucleic acid:

- (a) synthesizing cDNA molecules for each allele of each gene locus to be sequenced, wherein said synthesis employs a locus-specific oligonucleotide primer that anneals to a conserved region of each allele of each said gene locus.

86. A method for rapid automated determination of the genotype at one or more polymorphic gene locus of a subject in a sample containing subject nucleic acid comprising:

- (a) isolating nucleic acid from said sample with an RNA/DNA extractor;
- (b) amplifying said nucleic acid by polymerase chain reaction using a thermocycler to generate a polymerase chain reaction product for each allele of each gene locus to be sequenced, all of said alleles for each gene locus and chromosome to be sequenced being amplified with at least one conserved oligonucleotide primer and at least one conserved oligonucleotide primer pair and at least one of said alleles for each gene locus and chromosome to be sequenced being amplified with at least one conserved oligonucleotide primer and one non-conserved oligonucleotide primer;
- (c) sequencing directly each polymerase chain reaction product for each allele at each gene locus of each chromosome in an automated sequencing apparatus with a sequencing enzyme and a conserved primer specific for each locus to be sequenced; and

98

- (d) analyzing each sequenced polymerase chain reaction product to determine the genotype of said subject with a computer having a data base with allelic sequence information to compare the sequence of each allele of each gene locus sequenced to known sequences for each such gene locus followed by comparing the sequence of each allele of each gene locus amplified with a conserved/non-conserved oligonucleotide primer pair to the nucleotide sequence of each allele of such gene locus amplified with a conserved oligonucleotide primer pair.




FIG. 1B




FIG. 1A

SUBSTITUTE SHEET

SUBSTITUTE SHEET

4/4





FIG. 1C

SUBSTITUTE SHEET

FIG. 2A

4/4



SUBSTITUTE SHEET

